

Ethylenediaminetetraacetate (EDTA) disodium dihydrate

Safety Data Sheet

E0511

according to Regulation (EU) 2015/830

Version: 2.0

Date of issue: 13/09/2011 Revision date: 14/04/2017

SECTION 1: Identification of the substance/mixture and of the company/undertaking

1.1. Product identifier

Product form : Substance
 Substance name : Ethylenediaminetetraacetate (EDTA) disodium dihydrate
 EC-No. : 205-358-3
 CAS-No. : 6381-92-6
 REACH registration No : 01-2119486775-20-XXXX
 Product code : E0511
 Formula : C₁₀H₁₄N₂O₂Na₂ · 2H₂O
 Synonyms : EDTA · Na₂ salt dihydrate / EDTA disodium salt / Edetate disodium salt dihydrate

1.2. Relevant identified uses of the substance or mixture and uses advised against

1.2.1. Relevant identified uses

Main use category : Professional use
 Industrial/Professional use spec : For professional use only. Duchefa Biochemie B.V. products are intended only for "in vitro laboratory" research purposes

1.2.2. Uses advised against

No additional information available

1.3. Details of the supplier of the safety data sheet

Supplier

Duchefa Biochemie B.V.
 A. Hofmanweg 71
 2031 BH Haarlem - The Netherlands
 T +31(0)23-5319093 - F +31(0)23-5318027
info@duchefa.nl

1.4. Emergency telephone number

Emergency number : Supplier contact information:
 +31(0)23-5319093 (M-F 09:00-17:00)
 +31(0)6-30109355 (outside office hours)

Organisation/Company	Address	Comment
World Health Organization world directory of poison centres	http://apps.who.int/poisoncentres/	Consult website for a local poison centre

SECTION 2: Hazards identification

2.1. Classification of the substance or mixture

Classification according to Regulation (EC) No. 1272/2008 [CLP]

Acute toxicity (inhal.), Category 4 H332
 Specific target organ toxicity — Repeated exposure, Category 2 H373

Full text of hazard classes and H-statements : see section 16

Ethylenediaminetetraacetate (EDTA) disodium dihydrate

E0511

Safety Data Sheet

according to Regulation (EU) 2015/830

Adverse physicochemical, human health and environmental effects

No additional information available

2.2. Label elements

Labelling according to Regulation (EC) No. 1272/2008 [CLP]

Hazard pictograms (CLP)



GHS07

GHS08

Signal word (CLP)

: Warning

Hazard statements (CLP)

: H332 - Harmful if inhaled
H373 - May cause damage to organs through prolonged or repeated exposure

Precautionary statements (CLP)

: P260 - Do not breathe dust, fume, gas, mist, spray, vapours
P304+P340 - IF INHALED: remove victim to fresh air and keep at rest in a position comfortable for breathing
P312 - Call a POISON CENTER or doctor/physician if you feel unwell

2.3. Other hazards

No additional information available

SECTION 3: Composition/information on ingredients

3.1. Substances

Name : Ethylenediaminetetraacetate (EDTA) disodium dihydrate

CAS-No. : 6381-92-6

EC-No. : 205-358-3

Name	Product identifier	%
Ethylenediaminetetraacetate disodium dihydrate	(CAS-No.) 6381-92-6 (EC-No.) 205-358-3	>= 99

Full text of H-statements: see section 16

3.2. Mixtures

Not applicable

SECTION 4: First aid measures

4.1. Description of first aid measures

First-aid measures after inhalation : Move to fresh air. Get medical advice/attention if you feel unwell.

First-aid measures after skin contact : Wash with plenty of soap and water.

First-aid measures after eye contact : Rinse with plenty of water.

First-aid measures after ingestion : Rinse mouth. Do not induce vomiting. Seek medical attention if ill effect or irritation develops.

4.2. Most important symptoms and effects, both acute and delayed

No additional information available

Ethylenediaminetetraacetate (EDTA) disodium dihydrate

E0511

Safety Data Sheet

according to Regulation (EU) 2015/830

4.3. Indication of any immediate medical attention and special treatment needed

No additional information available

SECTION 5: Fire-fighting measures

5.1. Extinguishing media

Suitable extinguishing media : Water spray, carbon dioxide (CO₂), dry chemical powder, foam.

Unsuitable extinguishing media : Do not use a heavy water stream.

5.2. Special hazards arising from the substance or mixture

Hazardous decomposition products in case of fire : Under fire conditions, hazardous fumes will be present. CO_x. NO_x.

5.3. Advice for firefighters

Precautionary measures fire : Wear proper protective equipment.

Firefighting instructions : Prevent fire fighting water from entering the environment.

SECTION 6: Accidental release measures

6.1. Personal precautions, protective equipment and emergency procedures

6.1.1. For non-emergency personnel

Measures in case of dust release : Use good housekeeping practices to avoid rendering dust airborne.

6.1.2. For emergency responders

No additional information available

6.2. Environmental precautions

Prevent entry to sewers and public waters.

6.3. Methods and material for containment and cleaning up

Methods for cleaning up : Sweep up dry powder and dispose properly.

6.4. Reference to other sections

No additional information available

SECTION 7: Handling and storage

7.1. Precautions for safe handling

Precautions for safe handling : Open only in well ventilated area.

7.2. Conditions for safe storage, including any incompatibilities

Storage conditions : Store in dry, well-ventilated area. Store at room temperature.

7.3. Specific end use(s)

For professional use only. Duchefa Biochemie B.V. products are intended only for "in vitro laboratory" research purposes.

SECTION 8: Exposure controls/personal protection

8.1. Control parameters

Ethylenediaminetetraacetate (EDTA) disodium dihydrate (6381-92-6)
--

DNEL/DMEL (Workers)

Ethylenediaminetetraacetate (EDTA) disodium dihydrate

E0511

Safety Data Sheet

according to Regulation (EU) 2015/830

Ethylenediaminetetraacetate (EDTA) disodium dihydrate (6381-92-6)	
Acute - local effects, inhalation	3 mg/m ³
Long-term - local effects, inhalation	1,5 mg/m ³
PNEC (Water)	
PNEC aqua (freshwater)	2,2 mg/l
PNEC aqua (marine water)	0,22 mg/l
PNEC aqua (intermittent, freshwater)	1,2 mg/l
PNEC (Soil)	
PNEC soil	0,72 mg/kg dwt
PNEC (STP)	
PNEC sewage treatment plant	43 mg/l

8.2. Exposure controls

Hand protection:

Type	Material	Permeation	Thickness (mm)	Standard
Gloves	Nitrile rubber (NBR)	6 (> 480 minutes)	0,11	EN 374

Eye protection:

Chemical goggles or safety glasses

Skin and body protection:

Wear suitable protective clothing

Respiratory protection:

Where excessive dust may result, wear approved mask. Type P1 (EN 143)

SECTION 9: Physical and chemical properties

9.1. Information on basic physical and chemical properties

Physical state	: Solid
Appearance	: Crystals. Powder.
Molecular mass	: 372,2 g/mol
Colour	: White.
Odour	: Characteristic.
Odour threshold	: No data available
pH	: 4 (≥ 5) (aqueous solution 1%, 20 °C)
Relative evaporation rate (butylacetate=1)	: No data available
Melting point	: 252 °C
Freezing point	: No data available
Boiling point	: No data available
Flash point	: No data available
Auto-ignition temperature	: > 400 °C (EU Method A.16)
Decomposition temperature	: > 200 °C
Flammability (solid, gas)	: No data available
Vapour pressure	: No data available
Relative vapour density at 20 °C	: No data available
Relative density	: No data available
Density	: 1,77 g/cm ³ at 20 °C

Ethylenediaminetetraacetate (EDTA) disodium dihydrate

E0511

Safety Data Sheet

according to Regulation (EU) 2015/830

Solubility : Water: 108 g/l at 20 °C

Log Pow : No data available

Viscosity, kinematic : No data available

Viscosity, dynamic : No data available

Explosive properties : No data available

Oxidising properties : No data available

Explosive limits : No data available

9.2. Other information

Bulk density : 618 kg/m³

SECTION 10: Stability and reactivity

10.1. Reactivity

Stable under normal conditions of storage, handling and use.

10.2. Chemical stability

Stable under normal conditions.

10.3. Possibility of hazardous reactions

Stable under normal conditions of storage, handling and use.

10.4. Conditions to avoid

Direct sunlight. Avoid dust formation. Take precautionary measures against static discharges.

10.5. Incompatible materials

Metals. Strong acids. Strong bases. Oxidising agents.

10.6. Hazardous decomposition products

Carbon oxides (CO, CO₂).

SECTION 11: Toxicological information

11.1. Information on toxicological effects

Acute toxicity : Inhalation: Harmful if inhaled.

Ethylenediaminetetraacetate (EDTA) disodium dihydrate (6381-92-6)	
LD50 oral rat	> 2000 mg/kg
Skin corrosion/irritation	: Not classified pH: 4 (≥ 5) (aqueous solution 1%, 20 °C)
Serious eye damage/irritation	: Not classified pH: 4 (≥ 5) (aqueous solution 1%, 20 °C)
Respiratory or skin sensitisation	: Not classified
Germ cell mutagenicity	: Not classified
Carcinogenicity	: Not classified
Reproductive toxicity	: Not classified
STOT-single exposure	: Not classified
STOT-repeated exposure	: May cause damage to organs through prolonged or repeated exposure.

Ethylenediaminetetraacetate (EDTA) disodium dihydrate

E0511

Safety Data Sheet

according to Regulation (EU) 2015/830

Aspiration hazard : Not classified

SECTION 12: Ecological information

12.1. Toxicity

Ethylenediaminetetraacetate (EDTA) disodium dihydrate (6381-92-6)	
LC50 fish 1	> 100 mg/l <i>Lepomis macrochirus</i> (Bluegill)
EC50 Daphnia 1	> 100 mg/l EC50 48h - <i>Daphnia magna</i> [mg/l]
EC50 72h algae (1)	> 100 mg/l (<i>scenedemus subspicatus</i>)

12.2. Persistence and degradability

No additional information available

12.3. Bioaccumulative potential

Ethylenediaminetetraacetate (EDTA) disodium dihydrate (6381-92-6)	
Bioconcentration factor (BCF REACH)	1,8

12.4. Mobility in soil

No additional information available

12.5. Results of PBT and vPvB assessment

No additional information available

12.6. Other adverse effects

No additional information available

SECTION 13: Disposal considerations

13.1. Waste treatment methods

Waste treatment methods : Dispose in a safe manner in accordance with local/national regulations.
Avoid release to the environment.

SECTION 14: Transport information

In accordance with ADR / IATA / IMDG

ADR	IMDG	IATA
14.1. UN number		
Not applicable	Not applicable	Not applicable
14.2. UN proper shipping name		
Not applicable	Not applicable	Not applicable
Not applicable	Not applicable	Not applicable
14.3. Transport hazard class(es)		
Not applicable	Not applicable	Not applicable
Not applicable	Not applicable	Not applicable
14.4. Packing group		
Not applicable	Not applicable	Not applicable
14.5. Environmental hazards		
Not applicable	Not applicable	Not applicable
No supplementary information available		

Ethylenediaminetetraacetate (EDTA) disodium dihydrate

E0511

Safety Data Sheet

according to Regulation (EU) 2015/830

14.6. Special precautions for user

- Overland transport

Not applicable

- Transport by sea

Not applicable

- Air transport

Not applicable

14.7. Transport in bulk according to Annex II of Marpol and the IBC Code

Not applicable

SECTION 15: Regulatory information

15.1. Safety, health and environmental regulations/legislation specific for the substance or mixture

15.1.1. EU-Regulations

No REACH Annex XVII restrictions

Ethylenediaminetetraacetate (EDTA) disodium dihydrate is not on the REACH Candidate List

Ethylenediaminetetraacetate (EDTA) disodium dihydrate is not on the REACH Annex XIV List

15.1.2. National regulations

Ensure all national/local regulations are observed

Germany

VwVwS Annex reference : Water hazard class (WGK) 2, hazard to waters (KBwS-Beschluss; ID No. 104)

12th Ordinance Implementing the Federal Immission Control Act - 12.BImSchV : Is not subject of the 12. BImSchV (Hazardous Incident Ordinance)

Netherlands

SZW-lijst van kankerverwekkende stoffen : The substance is not listed

SZW-lijst van mutagene stoffen : The substance is not listed

NIET-limitatieve lijst van voor de voortplanting giftige stoffen – Borstvoeding : The substance is not listed

NIET-limitatieve lijst van voor de voortplanting giftige stoffen – Vruchtbaarheid : The substance is not listed

NIET-limitatieve lijst van voor de voortplanting giftige stoffen – Ontwikkeling : The substance is not listed

Denmark

Recommendations Danish Regulation : Young people below the age of 18 years are not allowed to use the product
Pregnant/breastfeeding women working with the product must not be in direct contact with the product

15.2. Chemical safety assessment

No additional information available

Ethylenediaminetetraacetate (EDTA) disodium dihydrate

E0511

Safety Data Sheet

according to Regulation (EU) 2015/830

SECTION 16: Other information

Indication of changes:

1.1	REACH registration No	Added	
2.2	Hazard pictograms (CLP)	Added	
2.2	Signal word (CLP)	Added	
2.2	Hazard statements (CLP)	Added	
2.2	Precautionary statements (CLP)	Added	
8.2	Hand protection	Modified	Specified material, thickness, et cetera of gloves
9.1	Auto-ignition temperature	Added	

Abbreviations and acronyms:

ADN	European Agreement concerning the International Carriage of Dangerous Goods by Inland Waterways
CLP	Classification Labelling Packaging Regulation; Regulation (EC) No 1272/2008
DPD	Dangerous Preparations Directive 1999/45/EC
DSD	Dangerous Substances Directive 67/548/EEC
LD50	Median lethal dose
REACH	Registration, Evaluation, Authorisation and Restriction of Chemicals Regulation (EC) No 1907/2006
SDS	Safety Data Sheet
ADR	European Agreement concerning the International Carriage of Dangerous Goods by Road
IATA	International Air Transport Association
IMDG	International Maritime Dangerous Goods
PBT	Persistent Bioaccumulative Toxic
EC50	Median effective concentration
LC50	Median lethal concentration
NOAEL	No-Observed Adverse Effect Level
LOAEL	Lowest Observed Adverse Effect Level
BCF	Bioconcentration factor

Data sources : Supplier. Manufacturer. REGULATION (EC) No 1272/2008 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL of 16 December 2008 on classification, labelling and packaging of substances and mixtures, amending and repealing Directives 67/548/EEC and 1999/45/EC, and amending Regulation (EC) No 1907/2006.

Full text of H- and EUH-statements:

Acute Tox. 4 (Inhalation)	Acute toxicity (inhal.), Category 4
STOT RE 2	Specific target organ toxicity — Repeated exposure, Category 2
H332	Harmful if inhaled
H373	May cause damage to organs through prolonged or repeated exposure

SDS Biochemicals version 2017

This information is based on our current knowledge and is intended to describe the product for the purposes of health, safety and environmental requirements only. It should not therefore be construed as guaranteeing any specific property of the product