

FICHE TECHNIQUE

F0123 : 5-FLUOROURACIL (5-FU)

CAS number 51-21-8

Molecular weight $C_4H_3FN_2O_2 = 130.1$

Assay > 98 %

Storage condition

- Store at 2-8 °C
- Soluble in water (10 g/l)
- Loss on drying < 0.5%

Hazard statements

- **H301** Toxic if swallowed
- **H312** Harmful in contact with skin
- **H332** Harmful if inhaled
- **H350** May cause cancer
- **H360** May damage fertility or the unborn child

Precautionary statements

- **P201** Obtain special instructions before use
- **P281** Use personal protective equipment as required
- **P301** IF SWALLOWED:
- **P308** IF exposed or concerned:
- **P310** Immediately call a POISON CENTER or doctor/physician
- **P313** Get medical advice/attention

F0123.0001/ 1 g

F0123.0005/ 5 g