

FICHE TECHNIQUE

F0176 : 5-FLUORO OROTIC ACID (5-FOA)

CAS number 703-95-7

Molecular weight $C_5H_3FN_2O_4 = 174.1$

Assay > 98 %

Storage condition

- Store at -25 to -15 °C
- Soluble in ethanol
- Slightly soluble in water

Hazard statements

- **H302** Harmful if swallowed
- **H315** Causes skin irritation
- **H319** Causes serious eye irritation

Precautionary statements

- **P280** Wear protective gloves/protective clothing/eye protection/face protection
- **P301** IF SWALLOWED:
- **P305** IF IN EYES:
- **P312** Call a POISON CENTER or doctor/physician if you feel **unwell**
- **P351** Rinse continuously with water for several minutes

F0176.1000/ 1g

F0176.5000/ 5g