

FICHE TECHNIQUE

F0512: FERROUS SULPHATE HEPTAHYDRATE

CAS number 7782-63-0

Molecular weight $\text{FeSO}_4 \cdot 7\text{H}_2\text{O}$ = 278.0

Assay > 98 %

Storage condition

- Store at room temperature
- Soluble in water (± 400 g/l / 20°C)

Hazard statements

- **H302** Harmful if swallowed
- **H315** Causes skin irritation
- **H319** Causes serious eye irritation

Precautionary statements

- **P280** Wear protective gloves/protective clothing/eye protection/face protection
- **P301** IF SWALLOWED:
- **P305** IF IN EYES:
- **P312** Call a POISON CENTER or doctor/physician if you feel unwell
- **P351** Rinse continuously with water for several minutes

F0512.0250/ 250 g

F0512.1000/ 1 kg

