

## FICHE TECHNIQUE



### **F0619: FOLINATE CALCIUM PENTAHYDRATE**

CAS number 1492-18-8

Molecular weight  $C_{20}H_{21}N_7O_7Ca \cdot 5H_2O = 601.5$

Assay > 97 %

Storage condition

- Store dry at 2-8 °C
- Soluble in water

#### **Hazard statements**

- **H315** Causes skin irritation
- **H317** May cause an allergic skin reaction
- **H319** Causes serious eye irritation
- **H334** May cause allergy or asthma symptoms or breathing difficulties if inhaled
- **H335** May cause respiratory irritation

#### **Precautionary statements**

- **P261** Avoid breathing dust/fume/gas/mist/vapours/spray
- **P280** Wear protective gloves/protective clothing/eye protection/face protection
- **P302** IF ON SKIN:
- **P305** IF IN EYES:
- **P351** Rinse continuously with water for several minutes
- **P352** Wash with soap and water

F0619.0001/ 1 g