

## FICHE TECHNIQUE



### **F0919 : FLURIDON**

CAS number 59756-60-4

Molecular weight C<sub>19</sub>H<sub>14</sub>F<sub>3</sub>NO = 329.3

Assay > 99 %

#### **Storage condition**

- Store at room temperature
- Slightly soluble in methanol

#### **Hazard statements**

- **H311** Toxic in contact with skin
- **H411** Toxic to aquatic life with long lasting effects

#### **Precautionary statements**

- **P273** Avoid release to the environment
- **P280** Wear protective gloves/protective clothing/eye protection/face protection
- **P312** Call a POISON CENTER or doctor/physician if you feel unwell

F0919.0250/ 250 mg