

## F3001

### FORCEPS 23 cm length stainless steel



#### Description :

This stainless steel forceps, with a length of 23 cm, has been especially developed for the handling of plantlets in Plant Tissue culture.

By means of the long thin extended legs, the distance between the hand, the -sterile plant-lets and culture vessel has been lengthened significantly, here-by drastically reducing the risks of contamination.

With its long thin legs it is easy maneuvering in long narrow culture tubes and due to the length, the bot-tom can be reached with-out contacting the sterile rim of the tubes.

Its light weight and the required low pressure by hand to close the forceps give it a fine ergo-no-metric performance without fatiguing the hand.

- **F3001.0001** : 1 pcs

## F3003

### FORCEPS 30 cm length stainless steel

#### Description :

This extended forceps allows a maximum distance between the hand and the plantlets minimizing the risks of contamination to nil.

The forceps is specially designed to work in combination with Ergonomic scalpel handle, S3310.

While handed the length of both tools is about equal providing a symmetric and ergonomically balanced work situation.

- **F3003.0001** : 1 pcs