

FICHE TECHNIQUE



G0124 : GENTAMYCIN SULPHATE

CAS number 1405-41-0
Potency > 560 units / mg

Storage condition

- Store at room temperature
- Soluble in water (100 g/l / 20 °C)

Hazard statements

- **H315** Causes skin irritation
- **H317** May cause an allergic skin reaction
- **H319** Causes serious eye irritation
- **H334** May cause allergy or asthma symptoms or breathing difficulties if inhaled

Precautionary statements

- **P280** Wear protective gloves/protective clothing/eye protection/face protection
- **P305** IF IN EYES:
- **P313** Get medical advice/attention
- **P337** If eye irritation persists:
- **P351** Rinse continuously with water for several minutes

G0124.0001/ 1 g
G0124.0005/ 5 g
G0124.0010/ 10 g
G0124.0025/ 25 g