

FICHE TECHNIQUE



G0175 : G-418 DISULPHATE

CAS number 108321-42-2

Molecular weight $C_{20}H_{40}N_{4}O_{10} \cdot 2H_2SO_4 = 692.7$

Assay :> 650 microg / mg

Storage condition

- Store at 2-8 °C
- Soluble in water

Hazard statements

- **H312** Harmful in contact with skin
- **H332** Harmful if inhaled
- **H360** May damage fertility or the unborn child

Precautionary statements

- **P280** Wear protective gloves/protective clothing/eye protection/face protection
- **P302** IF ON SKIN:
- **P304** IF INHALED:
- **P308** IF exposed or concerned:
- **P313** Get medical advice/attention
- **P341** If breathing is difficult
- **P352** Wash with soap and water

G0175.0001/ 1 g

G0175.0005/ 5 g