

FICHE TECHNIQUE

G0404: GRESSHOF & DOY (DBM2) VITAMIN MIXTURE

Vitamins

mg/l

µM

Glycine

4.00

53.28

myo-Inositol

100.00

554.94

Nicotinic acid

1.00

8.12

Pyridoxine HCl

1.00

4.86

Thiamine HCl

10.00

29.65

Gresshoff P.M. et al., Haploid Arabidopsis thaliana callus and plants from anther culture. Aust. J. Biol. Sci, 25, 259 (1972). Gresshoff P.M. et al. , Derivation of a haploid cell line from Vitis vinifera and the importance of the stage of meiotic development of anthers for haploid culture of this and other genera, Z. Pflanzenphysiol. 73, 132-141, (1974).

G0404.0100/ >100 ml