

FICHE TECHNIQUE

G0415: GAMBORG B5 VITAMIN MIXTURE

Vitamins

mg/l

µM

myo-Inositol

100.00

554.94

Nicotinic acid

1.00

8.12

Pyridoxine HCl

1.00

4.86

Thiamine HCl

10.00

29.65

Gamborg O.L., Miller R.A., Ojima K., Nutrient requirement of suspensions cultures of soybean root cells. Exp. Cell Res., 50, 151 (1968).

G0415.0100/ >100 ml

G0415.0250/ >250 ml