

# Gibberelic acid 3

## Safety Data Sheet

**G0907**

according to Regulation (EU) 2015/830

Version: 2.0

Date of issue: 07/11/2011    Revision date 06/07/2018  
Supersedes 06/03/2012

### SECTION 1: Identification of the substance/mixture and of the company/undertaking

#### 1.1. Product identifier

Product form : Substance  
Trade name : Gibberelic acid 3  
EC-No. : 201-001-0  
CAS-No. : 77-06-5  
Product code : G0907  
Formula : C<sub>19</sub>H<sub>22</sub>O<sub>6</sub>  
Synonyms : GA3  
Product group : Raw material

#### 1.2. Relevant identified uses of the substance or mixture and uses advised against

##### 1.2.1. Relevant identified uses

Main use category : Professional use  
Industrial/Professional use spec : For professional use only. Duchefa Biochemie B.V. products are intended only for "in vitro laboratory" research purposes.

##### 1.2.2. Uses advised against

No additional information available

#### 1.3. Details of the supplier of the safety data sheet

##### Supplier

Duchefa Biochemie B.V.  
A. Hofmanweg 71  
2031 BH Haarlem - The Netherlands  
T +31(0)23-5319093 - F +31(0)23-5318027  
[info@duchefa.nl](mailto:info@duchefa.nl)

#### 1.4. Emergency telephone number

Emergency number : Supplier contact information:  
+31(0)23-5319093 (M-F 09:00-17:00)  
+31(0)6-30109355 (outside office hours)

Organisation/Company	Address	Comment
World Health Organization world directory of poison centres	<a href="http://apps.who.int/poisoncentres/">http://apps.who.int/poisoncentres/</a>	Consult website for a local poison centre

### SECTION 2: Hazards identification

#### 2.1. Classification of the substance or mixture

##### Classification according to Regulation (EC) No. 1272/2008 [CLP]

Not classified

##### Adverse physicochemical, human health and environmental effects

No additional information available

## Safety Data Sheet

according to Regulation (EU) 2015/830

### 2.2. Label elements

#### Labelling according to Regulation (EC) No. 1272/2008 [CLP]

No labelling applicable

### 2.3. Other hazards

No additional information available

## SECTION 3: Composition/information on ingredients

### 3.1. Substances

Name : Gibberelic acid 3  
CAS-No. : 77-06-5  
EC-No. : 201-001-0

Name	Product identifier	%
Gibberellic acid 3	(CAS-No.) 77-06-5 (EC-No.) 201-001-0	>= 90

Full text of H-statements: see section 16

### 3.2. Mixtures

Not applicable

## SECTION 4: First aid measures

### 4.1. Description of first aid measures

First-aid measures after inhalation : Assure fresh air breathing.  
First-aid measures after skin contact : Wash skin with mild soap and water  
Seek medical attention if ill effect or irritation develops.  
First-aid measures after eye contact : Rinse with water  
Seek medical attention if ill effect or irritation develops.  
First-aid measures after ingestion : Rinse mouth  
Seek medical attention if ill effect or irritation develops.

### 4.2. Most important symptoms and effects, both acute and delayed

No additional information available

### 4.3. Indication of any immediate medical attention and special treatment needed

No additional information available

## SECTION 5: Firefighting measures

### 5.1. Extinguishing media

Suitable extinguishing media : Alcohol resistant foam.  
Dry chemical powder.  
Carbon dioxide (CO<sub>2</sub>).  
Water spray.

### 5.2. Special hazards arising from the substance or mixture

Hazardous decomposition products in case of fire : When exposed to heat, may decompose liberating hazardous gases  
- CO<sub>x</sub>.

### 5.3. Advice for firefighters

- Firefighting instructions : Use water spray or fog for cooling exposed containers.  
Protection during firefighting : Wear recommended personal protective equipment.

## SECTION 6: Accidental release measures

### 6.1. Personal precautions, protective equipment and emergency procedures

#### 6.1.1. For non-emergency personnel

Measures in case of dust release : Avoid dust formation.

#### 6.1.2. For emergency responders

No additional information available

### 6.2. Environmental precautions

Prevent entry to sewers and public waters.

### 6.3. Methods and material for containment and cleaning up

Methods for cleaning up : Sweep up dry powder and dispose properly.

### 6.4. Reference to other sections

No additional information available

## SECTION 7: Handling and storage

### 7.1. Precautions for safe handling

Precautions for safe handling : Handle in accordance with good industrial hygiene and safety practice.  
Avoid dust formation.

### 7.2. Conditions for safe storage, including any incompatibilities

Storage conditions : Store in dry, well-ventilated area  
Powder: store at room temperature  
Liquid: 2 - 8 °C  
Keep container closed when not in use.

### 7.3. Specific end use(s)

For professional use only. Duchefa Biochemie B.V. products are intended only for "in vitro laboratory" research purposes.

## SECTION 8: Exposure controls/personal protection

### 8.1. Control parameters

No additional information available

### 8.2. Exposure controls

#### Hand protection:

Type	Material	Permeation	Thickness (mm)	Standard
Gloves	Nitrile rubber (NBR)	6 (> 480 minutes)	0,11	EN 374

#### Eye protection:

Safety glasses (to European standard EN 166 or equivalent)

#### Skin and body protection:

Wear suitable protective clothing

**Respiratory protection:**

Where excessive dust may result, wear approved mask. Type P1 (EN 143)

**SECTION 9: Physical and chemical properties****9.1. Information on basic physical and chemical properties**

Physical state	: Solid
Appearance	: Crystalline powder.
Molecular mass	: 346,4 g/mol
Colour	: White.
Odour	: Mild. Odourless.
Odour threshold	: No data available
pH	: No data available
Relative evaporation rate (butylacetate=1)	: No data available
Melting point	: 232 - 235 °C
Freezing point	: No data available
Boiling point	: No data available
Flash point	: No data available
Auto-ignition temperature	: No data available
Decomposition temperature	: No data available
Flammability (solid, gas)	: No data available
Vapour pressure	: No data available
Relative vapour density at 20 °C	: No data available
Relative density	: 0,75 - 0,85
Solubility	: Soluble in ethanol.
Log Pow	: No data available
Viscosity, kinematic	: No data available
Viscosity, dynamic	: No data available
Explosive properties	: No data available
Oxidising properties	: No data available
Explosive limits	: No data available

**9.2. Other information**

No additional information available

**SECTION 10: Stability and reactivity****10.1. Reactivity**

Stable under normal conditions.

**10.2. Chemical stability**

Stable under normal conditions of storage, handling and use.

**10.3. Possibility of hazardous reactions**

No additional information available

**10.4. Conditions to avoid**

Moisture.

### 10.5. Incompatible materials

Strong oxidizers.

### 10.6. Hazardous decomposition products

According to process conditions, hazardous decomposition products may be generated.

- CO<sub>x</sub>.

## SECTION 11: Toxicological information

### 11.1. Information on toxicological effects

Acute toxicity : Not classified

Gibberelic acid 3 (77-06-5)	
LD50 oral rat	6300 mg/kg
LD50 dermal rabbit	> 2000 mg/kg

Skin corrosion/irritation : Not classified

Serious eye damage/irritation : Not classified

Respiratory or skin sensitisation : Not classified

Germ cell mutagenicity : Not classified

Carcinogenicity : Not classified

Reproductive toxicity : Not classified

STOT-single exposure : Not classified

STOT-repeated exposure : Not classified

Aspiration hazard : Not classified

## SECTION 12: Ecological information

### 12.1. Toxicity

Gibberelic acid 3 (77-06-5)	
LC50 fish 1	> 150 ppm <i>Oncorhynchus mykiss</i> (Rainbow trout)

### 12.2. Persistence and degradability

Gibberelic acid 3 (77-06-5)	
BOD (% of ThOD)	≈ 99 % ThOD

### 12.3. Bioaccumulative potential

No additional information available

### 12.4. Mobility in soil

No additional information available

### 12.5. Results of PBT and vPvB assessment

No additional information available

### 12.6. Other adverse effects

Additional information : Prevent entry to sewers and public waters

**SECTION 13: Disposal considerations****13.1. Waste treatment methods**

Waste treatment methods : Dispose of this material and its container at hazardous or special waste collection point.

**SECTION 14: Transport information**

In accordance with ADR / IATA / IMDG

ADR	IMDG	IATA
<b>14.1. UN number</b>		
Not regulated	Not regulated	Not regulated
<b>14.2. UN proper shipping name</b>		
Not regulated	Not regulated	Not regulated
Not regulated	Not regulated	Not regulated
<b>14.3. Transport hazard class(es)</b>		
Not regulated	Not regulated	Not regulated
Not regulated	Not regulated	Not regulated
<b>14.4. Packing group</b>		
Not regulated	Not regulated	Not regulated
<b>14.5. Environmental hazards</b>		
Not regulated	Not regulated	Not regulated
No supplementary information available		

**14.6. Special precautions for user****. Overland transport**

Not regulated

**- Transport by sea**

Not regulated

**- Air transport**

Not regulated

**14.7. Transport in bulk according to Annex II of Marpol and the IBC Code**

Not applicable

**SECTION 15: Regulatory information****15.1. Safety, health and environmental regulations/legislation specific for the substance or mixture****15.1.1. EU-Regulations**

No REACH Annex XVII restrictions

Gibberelic acid 3 is not on the REACH Candidate List

Gibberelic acid 3 is not on the REACH Annex XIV List

**15.1.2. National regulations**

Ensure all national/local regulations are observed.

**Germany**

12th Ordinance Implementing the Federal Immission Control Act -  
12.BImSchV

: Is not subject of the 12. BImSchV  
(Hazardous Incident Ordinance)

## Safety Data Sheet

according to Regulation (EU) 2015/830

**Netherlands**

SZW-lijst van kankerverwekkende stoffen	: The substance is not listed
SZW-lijst van mutagene stoffen	: The substance is not listed
NIET-limitatieve lijst van voor de voortplanting giftige stoffen – Borstvoeding	: The substance is not listed
NIET-limitatieve lijst van voor de voortplanting giftige stoffen – Vruchtbaarheid	: The substance is not listed
NIET-limitatieve lijst van voor de voortplanting giftige stoffen – Ontwikkeling	: The substance is not listed

**15.2. Chemical safety assessment**

No additional information available

**SECTION 16: Other information**

Indication of changes:

1.4	Emergency number	Modified	
8.2	Hand protection	Modified	Specified material, thickness, et cetera of gloves

Abbreviations and acronyms:

ATE	Acute Toxicity Estimate
ADR	European Agreement concerning the International Carriage of Dangerous Goods by Road
BCF	Bioconcentration factor
CLP	Classification Labelling Packaging Regulation; Regulation (EC) No 1272/2008
DPD	Dangerous Preparations Directive 1999/45/EC
DSD	Dangerous Substances Directive 67/548/EEC
IATA	International Air Transport Association
IMDG	International Maritime Dangerous Goods
LC50	Median lethal concentration
LD50	Median lethal dose
LOAEL	Lowest Observed Adverse Effect Level
NOAEC	No-Observed Adverse Effect Concentration
PBT	Persistent Bioaccumulative Toxic
REACH	Registration, Evaluation, Authorisation and Restriction of Chemicals Regulation (EC) No 1907/2006
SDS	Safety Data Sheet

Data sources : ECHA (European Chemicals Agency). Supplier.

Full text of H- and EUH-statements:

Eye Irrit. 2	Serious eye damage/eye irritation, Category 2
H319	Causes serious eye irritation.

SDS Biochemicals version 2018

*This information is based on our current knowledge and is intended to describe the product for the purposes of health, safety and environmental requirements only. It should not therefore be construed as guaranteeing any specific property of the product*