

FICHE TECHNIQUE



G1375: GUANIDINE HCL

CAS number 50-01-1

Molecular weight $\text{CH}_5\text{N}_3\text{HCl}$ = 95.5

Assay > 99.5 %

Melting point 185 - 188

Moisture content < 0.2 %

Storage condition

- Store at room temperature
- Soluble in water (2280 g/l / 30 °C)
- Hygroscopic

Hazard statements

- **H302** Harmful if swallowed
- **H315** Causes skin irritation
- **H319** Causes serious eye irritation

Precautionary statements

- **P280** Wear protective gloves/protective clothing/eye protection/face protection
- **P302** IF ON SKIN:
- **P305** IF IN EYES:
- **P351** Rinse continuously with water for several minutes
- **P352** Wash with soap and water

G1375.0100/ 100 g

G1375.0250/ 250 g

G1375.0500/ 500 g

G1375.1000/ 1 kg