

## G3301

### Glass Bead Sterilizer 220-230V. / EU plug + 1x filling of Glass beads

Weight 3.5 kg  
Outside dimensions 15 x 12 x 15 cm  
Tube dimensions (diameter x height) 5.5 x 12 cm  
Operating temperature 275 C.  
Thermostat controlled 200 W / 220 Volt, 110 Volt upon request  
Including glass beads

#### Technical info

Storage condition

Outside dimensions: 15x12x15 cm, Weight: 3.5 kg  
Operating temperature 275 °C  
200 W/220 Volt or 110 Volt

- **G3301.0001** : 1 pcs

## G3302

### Glass Beads for Sterilizer G3301 (jar for 1 filling 500g)

Diameter 1.5 - 2 mm

#### Technical info

Storage condition

Diameter: 1.5 - 2 mm  
to be used in a glass bead sterilizer

- **G3302.0500** : 500 g

