

FICHE TECHNIQUE

H0168 : 8-HYDROXYQUINOLINE

CAS number 148-24-3

Molecular weight C₉H₇NO = 145.2

Assay > 99 %

Storage condition

- Store dark at room temperature
- Insoluble in water

Hazard statements

- **H302** Harmful if swallowed
- **H332** Harmful if inhaled

Precautionary statements

- **P281** Use personal protective equipment as required
- **P301** IF SWALLOWED:
- **P304** IF INHALED:
- **P312** Call a POISON CENTER or doctor/physician if you feel unwell
- **P340** Remove victim to fresh air and keep at rest in a position comfortable for breathing

H0168.0025/ 25g

H0168.0100/ 100 g