

FICHE TECHNIQUE



H0213: HELLER MEDIUM

Micro Elements

mg/l

μM

AlCl₃.6H₂O

0.054

0.22

CuSO₄.5H₂O

0.03

0.12

FeCl₃.6H₂O

1.00

3.70

H₃BO₃

6.20

100.27

KI

0.015

0.09

MnSO₄.H₂O

0.08

0.47

NiCl₂.6H₂O

0.025

0.14

ZnSO₄.7H₂O

1.00

3.48

Macro Elements

mg/l

mM

CaCl₂

56.62

0.51

KCl

750.00

10.06

MgSO₄

121.56

1.01

NaNO₃

600.00

7.06

NaH₂PO₄

108.70

0.91

Heller R., Ann. Sci. Nat. Bot. Biol. Veg. 11th Ser., 14, 1 (1953).

Hazard statements

- **H315** Causes skin irritation
- **H319** Causes serious eye irritation
- **H335** May cause respiratory irritation

Precautionary statements

- **P261** Avoid breathing dust/fume/gas/mist/vapours/spray
- **P302** IF ON SKIN:
- **P305** IF IN EYES:
- **P351** Rinse continuously with water for several minutes
- **P352** Wash with soap and water