

FICHE TECHNIQUE



I0902 : INDOLE-3-BUTYRIC ACID (IBA)

CAS number 133-32-4
Molecular weight $C_{12}H_{13}NO_2 = 203.2$
Assay > 98.0 %
Melting point 122 - 124 C.

Storage condition

- Store at 2-8 °C
- Soluble in ethanol and 1N NaOH

Hazard statements

- **H301** Toxic if swallowed
- **H315** Causes skin irritation
- **H319** Causes serious eye irritation
- **H335** May cause respiratory irritation

Precautionary statements

- **P280** Wear protective gloves/protective clothing/eye protection/face protection
- **P301** IF SWALLOWED:
- **P305** IF IN EYES:
- **P310** Immediately call a POISON CENTER or doctor/physician
- **P351** Rinse continuously with water for several minutes

I0902.0005/ 5 g
I0902.0025/ 25 g