

FICHE TECHNIQUE



K0215: KNUDSON C ORCHID MEDIUM

Micro Elements

mg/l

μM

FeSO₄.7H₂O

25.00

89.92

MnSO₄.H₂O

5.68

33.61

Macro Elements

mg/l

mM

Ca(NO₃)₂

241.30

1.43

KCl

250.00

3.35

KH₂PO₄

250.00

1.84

MgSO₄

122.15

1.02

NH₄NO₃

500.00

6.25

(NH₄)₂SO₄

500.00

3.78

Morel, G.M., Cymb. Soc. News, 20, (1965).

Hazard statements

- **H319** Causes serious eye irritation

Precautionary statements

- **P280** Wear protective gloves/protective clothing/eye protection/face protection
- **P305** IF IN EYES:
- **P351** Rinse continuously with water for several minutes

K0215.0001/ 1 l

K0215.0005/ 5 l

K0215.0010/ 10 l

K0215.0025/ 25 l

K0215.0050/ 50 l