

Safety Data Sheet

according to Regulation (EU) 2015/830 Date of issue: 30/03/2012 Revision date 17/10/2018 K5120

Version: 2.0

# SECTION 1: Identification of the substance/mixture and of the company/undertaking

Product form	: Mixture
Trade name	: KBBC Medium
Product code	: K5120
Product group	: Blend

#### 1.2. Relevant identified uses of the substance or mixture and uses advised against

#### 1.2.1. Relevant identified uses

 Main use category
 : Professional use

 Industrial/Professional use spec
 : For professional use only. Duchefa Biochemie B.V. products are intended only for "in vitro laboratory" research purposes.

#### 1.2.2. Uses advised against

No additional information available

#### 1.3. Details of the supplier of the safety data sheet

#### Manufacturer

Duchefa Biochemie B.V. A. Hofmanweg 71 2031 BH Haarlem - The Netherlands T +31(0)23-5319093 - F +31(0)23-5318027 info@duchefa.nl

#### 1.4. Emergency telephone number

Emergency number

- : Supplier contact information:
  - +31(0)23-5319093 (M-F 09:00-17:00)
  - +31(0)6-30109355 (outside office hours)

Organisation/Company	Address	Comment
World Health Organization world directory of poison centres	http://apps.who.int/poisoncentres/	Consult website for a local poison centre

### **SECTION 2: Hazards identification**

#### 2.1. Classification of the substance or mixture

#### Classification according to Regulation (EC) No. 1272/2008 [CLP] Not classified

#### Adverse physicochemical, human health and environmental effects

No additional information available

#### 2.2. Label elements

#### Labelling according to Regulation (EC) No. 1272/2008 [CLP]

No labelling applicable

#### 2.3. Other hazards

Safety Data Sheet

according to Regulation (EU) 2015/830

# SECTION 3: Composition/information on ingredients

#### 3.1. Substances

Not applicable

#### 3.2. Mixtures

Name	Product identifier	%	Classification according to Regulation (EC) No. 1272/2008 [CLP]
Proteose Peptone	(CAS-No.) 73049-73-7	51,65	Not classified
Agar Agar	(CAS-No.) 9002-18-0 (EC-No.) 232-658-1	38,73	Not classified
Boric acid substance listed as REACH Candidate	(CAS-No.) 10043-35-3 (EC-No.) 233-139-2 (EC Index-No.) 005-007-00-2	3,87	Repr. 1B, H360FD
Dipotassium hydrogenphosphate	(CAS-No.) 7758-11-4 (EC-No.) 231-834-5	3,87	Not classified
Magnesium sulfate anhydrous	(CAS-No.) 7487-88-9 (EC-No.) 231-298-2	1,88	Not classified

#### **Specific concentration limits:**

Name	Product identifier	Specific concentration limits
Boric acid	(CAS-No.) 10043-35-3 (EC-No.) 233-139-2 (EC Index-No.) 005-007-00-2	(C >= 5,5) Repr. 1B, H360FD

Full text of H-statements: see section 16

### SECTION 4: First aid measures

#### 4.1. Description of first aid measures

First-aid measures general	: Seek medical attention if ill effect develops.
First-aid measures after inhalation	: Remove victim to fresh air.
First-aid measures after skin contact	: Wash skin with mild soap and water.
First-aid measures after eye contact	: Rinse with water.
First-aid measures after ingestion	: Rinse mouth.
4.2 Most important symptoms a	nd offects both acute and delayed

#### 4.2. Most important symptoms and effects, both acute and delayed

Symptoms/effects after skin contact	: Causes skin irritation.
Symptoms/effects after eye contact	: Redness, pain.

# 4.3. Indication of any immediate medical attention and special treatment needed

No additional information available

Safety Data Sheet

according to Regulation (EU) 2015/830

# **SECTION 5: Firefighting measures**

SECTION 5: Firefighting mea	asures
<b>5.1. Extinguishing media</b> Suitable extinguishing media	: Alcohol resistant foam. Dry chemical powder. Carbon dioxide (CO2). Water spray.
<b>5.2.</b> Special hazards arising from Hazardous decomposition products in case of fire	<ul> <li>the substance or mixture</li> <li>Under fire conditions, hazardous fumes will be present:</li> <li>- COx</li> <li>- NOx</li> <li>- SOx.</li> </ul>
<b>5.3.</b> Advice for firefighters Firefighting instructions Protection during firefighting	<ul><li>Prevent fire fighting water from entering the environment.</li><li>Wear proper protective equipment.</li></ul>
<b>SECTION 6: Accidental relea</b>	ase measures
6.1. Personal precautions, protect General measures	tive equipment and emergency procedures : Avoid raising powdered materials into airborne dust.
6.1.1. For non-emergency personne Emergency procedures	el : Wear suitable protective clothing.
<b>6.1.2. For emergency responders</b> No additional information available	
<b>6.2. Environmental precautions</b> Prevent entry to sewers and public wat	ers.
6.3. Methods and material for cor Methods for cleaning up	ntainment and cleaning up : Sweep up dry powder and dispose properly.
6.4. Reference to other sections No additional information available	
SECTION 7: Handling and st	torage
7.1. Precautions for safe handling	

Precautions for safe handling : Avoid dust formation. Handle in accordance with good industrial hygiene and safety procedures.

### 7.2. Conditions for safe storage, including any incompatibilities

Storage conditions

: Store at room temperature Store in dry, well-ventilated area Hygroscopic Keep container tightly closed and dry.

#### 7.3. Specific end use(s)

For professional use only. Duchefa Biochemie B.V. products are intended only for "in vitro laboratory" research purposes.

Safety Data Sheet

according to Regulation (EU) 2015/830

# **SECTION 8: Exposure controls/personal protection**

#### 8.1. Control parameters

No additional information available

#### 8.2. Exposure controls

#### Hand protection:

Туре	Material	Permeation	Thickness (mm)	Standard
Gloves	Nitrile rubber (NBR)	6 (> 480 minutes)	0,11	EN 374

#### Eye protection:

Safety glasses (to European standard EN 166 or equivalent)

#### Skin and body protection:

In case of possible repeated skin contact wear protective clothing

#### **Respiratory protection:**

Where excessive dust may result, wear approved mask. Type P1 (EN 143)

### **SECTION 9: Physical and chemical properties**

#### 9.1. Information on basic physical and chemical properties

Physical state Appearance Colour	<ul><li>Solid</li><li>Powder.</li><li>White to slightly yellow.</li></ul>
Odour Odour threshold pH Relative evaporation rate (butylacetate=1)	<ul> <li>Characteristic. Weak.</li> <li>No data available</li> <li>No data available</li> <li>No data available</li> </ul>
Melting point Freezing point Boiling point Flash point Auto-ignition temperature Decomposition temperature Flammability (solid, gas) Vapour pressure Relative vapour density at 20 °C Relative density Solubility	<ul> <li>No data available</li> <li>Readily soluble in water.</li> </ul>
Log Pow Viscosity, kinematic Viscosity, dynamic Explosive properties Oxidising properties Explosive limits	<ul> <li>No data available</li> </ul>

Safety Data Sheet

according to Regulation (EU) 2015/830

#### 9.2. Other information

No additional information available

### SECTION 10: Stability and reactivity

#### 10.1. Reactivity

Stable under normal conditions of storage, handling and use.

#### 10.2. Chemical stability

Stable under normal conditions.

#### 10.3. Possibility of hazardous reactions

No additional information available

#### **10.4.** Conditions to avoid

Moisture.

#### 10.5. Incompatible materials

Strong oxidizers.

#### 10.6. Hazardous decomposition products

Thermal decomposition generates :

- COx
- NOx
- SOx
- Phosphorus oxides.

## **SECTION 11: Toxicological information**

#### 11.1. Information on toxicological effects

Acute toxicity	: Not classified
Skin corrosion/irritation	: Not classified
Serious eye damage/irritation	: Not classified
Respiratory or skin sensitisation	: Not classified
Germ cell mutagenicity	: Not classified
Carcinogenicity	: Not classified
Reproductive toxicity	: Not classified
STOT-single exposure	: Not classified
STOT-repeated exposure	: Not classified
Aspiration hazard	: Not classified

# SECTION 12: Ecological information

#### 12.1. Toxicity

No additional information available

#### 12.2. Persistence and degradability

No additional information available

5/8

Safety Data Sheet

according to Regulation (EU) 2015/830

#### 12.3. Bioaccumulative potential

No additional information available

#### 12.4. Mobility in soil

No additional information available

#### 12.5. Results of PBT and vPvB assessment

KBBC Medium	
Component	
Boric acid (10043-35-3)	This substance/mixture does not meet the PBT criteria of REACH regulation, annex XIII This substance/mixture does not meet the vPvB criteria of REACH regulation, annex XIII

#### 12.6. Other adverse effects

# Additional information

: Prevent entry to sewers and public waters. Avoid release to the environment

### SECTION 13: Disposal considerations

#### 13.1. Waste treatment methods

- Waste treatment methods
- : Dispose in a safe manner in accordance with local/national regulations. Avoid release to the environment.

# **SECTION 14: Transport information**

In accordance with ADR / IATA / IMDG

ADR	IMDG	ΙΑΤΑ	
14.1. UN number			
Not regulated	Not regulated	Not regulated	
14.2. UN proper shipping name			
Not regulated	Not regulated	Not regulated	
Not regulated	Not regulated	Not regulated	
14.3. Transport hazard class(es)	14.3. Transport hazard class(es)		
Not regulated	Not regulated	Not regulated	
Not regulated	Not regulated	Not regulated	
14.4. Packing group			
Not regulated	Not regulated	Not regulated	
14.5. Environmental hazards			
Not regulated	Not regulated	Not regulated	
No supplementary information available			

#### 14.6. Special precautions for user

#### Overland transport

Not regulated

#### - Transport by sea

Not regulated

#### - Air transport

Not regulated

Safety Data Sheet

according to Regulation (EU) 2015/830

#### 14.7. Transport in bulk according to Annex II of Marpol and the IBC Code

Not applicable

### SECTION 15: Regulatory information

# 15.1. Safety, health and environmental regulations/legislation specific for the substance or mixture

#### 15.1.1. EU-Regulations

Contains no REACH substances with Annex XVII restrictions Contains a substance on the REACH candidate list in concentration  $\geq$  0.1% or with a lower specific limit: Boric acid (EC 233-139-2, CAS 10043-35-3)

Contains no REACH Annex XIV substances

#### 15.1.2. National regulations

Ensure all national/local regulations are observed.

#### Germany

Reference to AwSV : Water hazard class (WGK) 3. severe hazard to water (Classification according to AwSV, Annex 1) 12th Ordinance Implementing the Federal Immission Control Act -: Is not subject of the 12. BImSchV 12.BImSchV (Hazardous Incident Ordinance) **Netherlands** SZW-lijst van kankerverwekkende stoffen : Agar Agar is listed SZW-lijst van mutagene stoffen : Agar Agar is listed NIET-limitatieve lijst van voor de voortplanting giftige stoffen - Borstvoeding : None of the components are listed NIET-limitatieve lijst van voor de voortplanting giftige stoffen -: Boric acid is listed Vruchtbaarheid NIET-limitatieve lijst van voor de voortplanting giftige stoffen - Ontwikkeling : Boric acid is listed Denmark **Recommendations Danish Regulation** : The requirements from the Danish

#### 15.2. Chemical safety assessment

No additional information available

### **SECTION 16: Other information**

Indication of changes:

8.2	Hand protection	Modified	Specified material, thickness, et	
			cetera of gloves	

Abbreviations and acronyms:

ATE	Acute Toxicity Estimate
ADR	European Agreement concerning the International Carriage of Dangerous Goods by Road
BCF	Bioconcentration factor
CLP	Classification Labelling Packaging Regulation; Regulation (EC) No 1272/2008
DPD	Dangerous Preparations Directive 1999/45/EC
DSD	Dangerous Substances Directive 67/548/EEC
IATA	International Air Transport Association

The requirements from the Danish Working Environment Authorities regarding work with carcinogens must be followed during use and disposal

### Safety Data Sheet

according to Regulation (EU) 2015/830

IMDG	International Maritime Dangerous Goods
LC50	Median lethal concentration
LD50	Median lethal dose
LOAEL	Lowest Observed Adverse Effect Level
NOAEC No-Observed Adverse Effect Concentration	
PBT Persistent Bioaccumulative Toxic	
REACH	Registration, Evaluation, Authorisation and Restriction of Chemicals Regulation (EC) No 1907/2006
SDS Safety Data Sheet	

#### Data sources

: REGULATION (EC) No 1272/2008 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL of 16 December 2008 on classification, labelling and packaging of substances and mixtures, amending and repealing Directives 67/548/EEC and 1999/45/EC, and amending Regulation (EC) No 1907/2006. TNO (Netherlands Organisation for Applied Scientific Research).

#### Full text of H- and EUH-statements:

Repr. 1B	Reproductive toxicity, Category 1B
H360FD	May damage fertility. May damage the unborn child.

#### SDS Biochemicals version 2018

This information is based on our current knowledge and is intended to describe the product for the purposes of health, safety and environmental requirements only. It should not therefore be construed as guaranteeing any specific property of the product