

## FICHE TECHNIQUE



### **L0216: LINDEMANN ORCHID MEDIUM**

#### **Micro Elements**

**mg/l**

**μM**

AlCl<sub>3</sub>.6H<sub>2</sub>O

0.56

2.32

CuSO<sub>4</sub>.5H<sub>2</sub>O

0.02

0.08

FeCitrate

4.40

H<sub>3</sub>BO<sub>3</sub>

1.01

16.34

KI

0.10

0.60

MnSO<sub>4</sub>.H<sub>2</sub>O

0.05

0.31

NiCl<sub>2</sub>.6H<sub>2</sub>O

0.03

0.13

ZnSO<sub>4</sub>.7H<sub>2</sub>O

0.57

1.98

## Macro Elements

mg/l

mM

Ca(NO<sub>3</sub>)<sub>2</sub>

347.20

2.12

KH<sub>2</sub>PO<sub>4</sub>

135.00

0.99

KCl

1050.00

14.08

MgSO<sub>4</sub>

58.98

0.49

(NH<sub>4</sub>)<sub>2</sub>SO<sub>4</sub>

1000.00

7.57

Lindemann E.G.P., Amercan. Orch. Bull 39, 1002 (1970).

### Hazard statements

- **H315** Causes skin irritation
- **H319** Causes serious eye irritation
- **H335** May cause respiratory irritation

### Precautionary statements

- **P261** Avoid breathing dust/fume/gas/mist/vapours/spray
- **P302** IF ON SKIN:
- **P305** IF IN EYES:
- **P351** Rinse continuously with water for several minutes
- **P352** Wash with soap and water

L0216.0001/ 1 l

L0216.0010/ 10 l

L0216.0050/ 50 l