

FICHE TECHNIQUE

L1204 : LOW MELTING AGAROSE PPC

CAS number 9012-36-6
Gel strength 1.5 % > 920g/cm²
Gelling temperature 24 - 30 C.
Melting temperature approximately 65 C.
Sulphate < 0.12 %
Elektroendosmosis 0.12
Moisture approximately 5 %

Storage condition

- Store at room temperature
- Soluble in water, gel strength: > 1000 g/cm²
- Melting temperature: < 65 °C

L1204.0100/ 100 g
L1204.0250/ 250 g