

FICHE TECHNIQUE

L1349: D-LUCIFERIN

CAS number 2591-17-5

Molecular weight $C_{11}H_8N_2O_3S_2 = 280.3$

Assay

D-Luciferin, HPLC, chemical purity

D-Luciferin, optical purity

Contains 0.05 % acetic acid as antistatic

Storage condition

- Store at -25 to -15 °C
- Protect from light and moisture
- Contains 0.05% acetic acid as antistatic

L1349.0100/ 100 mg

L1349.0250/ 250 mg

L1349.0500/ 500 mg

L1349.1000/ 1 g