

Safety Data Sheet

L1704

Version: 2.0

according to Regulation (EU) 2015/830 Date of issue: 02/01/2013 Revision date 17/10/2018

SECTION 1: Identification of the substance/mixture and of the company/undertaking	
1.1.	Product identifier

Product form	: Mixture
Trade name	: LB Broth High salt
Product code	: L1704
Product group	: Blend

#### 1.2. Relevant identified uses of the substance or mixture and uses advised against

#### 1.2.1. Relevant identified uses

 Main use category
 : Professional use

 Industrial/Professional use spec
 : For professional use only. Duchefa Biochemie B.V. products are intended only for "in vitro laboratory" research purposes.

#### 1.2.2. Uses advised against

No additional information available

#### 1.3. Details of the supplier of the safety data sheet

#### Manufacturer

Duchefa Biochemie B.V. A. Hofmanweg 71 2031 BH Haarlem - The Netherlands T +31(0)23-5319093 - F +31(0)23-5318027 info@duchefa.nl

#### 1.4. Emergency telephone number

Emergency number

- : Supplier contact information:
  - +31(0)23-5319093 (M-F 09:00-17:00)
  - +31(0)6-30109355 (outside office hours)

Organisation/Company	Address	Comment
World Health Organization world directory of poison centres	http://apps.who.int/poisoncentres/	Consult website for a local poison centre

#### **SECTION 2: Hazards identification**

#### 2.1. Classification of the substance or mixture

#### Classification according to Regulation (EC) No. 1272/2008 [CLP] Not classified

#### Adverse physicochemical, human health and environmental effects

No additional information available

#### 2.2. Label elements

#### Labelling according to Regulation (EC) No. 1272/2008 [CLP]

No labelling applicable

#### 2.3. Other hazards

Safety Data Sheet

according to Regulation (EU) 2015/830

#### SECTION 3: Composition/information on ingredients

#### 3.1. Substances

Not applicable

#### 3.2. Mixtures

Name	Product identifier	%	Classification according to Regulation (EC) No. 1272/2008 [CLP]
Sodium chloride	(CAS-No.) 7647-14-5 (EC-No.) 231-598-3	40	Not classified
Tryptone	(CAS-No.) 73049-73-7	40	Not classified
Yeast extract	(CAS-No.) 8013-01-2 (EC-No.) 232-387-9	20	Not classified

#### SECTION 4: First aid measures

#### 4.1. Description of first aid measures

First-aid measures general	: Seek medical attention if ill effect develops.	
First-aid measures after inhalation	: Remove victim to fresh air.	
First-aid measures after skin contact	: Wash skin with mild soap and water.	
First-aid measures after eye contact	: Rinse with water.	
First-aid measures after ingestion : Rinse mouth.		
4.2. Most important symptoms and effects, both acute and delayed		

Symptoms/effects after skin contact	: Causes skin irritation.
Symptoms/effects after eye contact	: Redness, pain.

## **4.3.** Indication of any immediate medical attention and special treatment needed No additional information available

## SECTION 5: Firefighting measures

#### 5.1. Extinguishing media Suitable extinguishing media : Alcohol resistant foam. Dry chemical powder. Carbon dioxide (CO2).

Water spray.

# 5.2. Special hazards arising from the substance or mixture Hazardous decomposition products in case of fire : Under fire conditions, hazardous fumes will be present: - COx - NOx

#### - SOx.

#### 5.3. Advice for firefighters

Firefighting instructions	: Prevent fire fighting water from entering the environment.
Protection during firefighting	: Wear proper protective equipment.

Safety Data Sheet

according to Regulation (EU) 2015/830

#### **SECTION 6: Accidental release measures**

#### 6.1. Personal precautions, protective equipment and emergency procedures General measures : Avoid raising powdered materials into airborne dust. 6.1.1. For non-emergency personnel Emergency procedures : Wear suitable protective clothing. 6.1.2. For emergency responders No additional information available 6.2. **Environmental precautions** Prevent entry to sewers and public waters. Methods and material for containment and cleaning up 6.3. : Sweep up dry powder and dispose properly. Methods for cleaning up 6.4. **Reference to other sections** No additional information available SECTION 7: Handling and storage 7.1. Precautions for safe handling Precautions for safe handling : Avoid dust formation. Handle in accordance with good industrial hygiene and safety procedures. 7.2. Conditions for safe storage, including any incompatibilities Storage conditions : Store at room temperature Store in dry, well-ventilated area Hygroscopic

#### 7.3. Specific end use(s)

For professional use only. Duchefa Biochemie B.V. products are intended only for "in vitro laboratory" research purposes.

Keep container tightly closed and dry.

#### SECTION 8: Exposure controls/personal protection

#### 8.1. Control parameters

No additional information available

#### 8.2. Exposure controls

#### Hand protection:

Туре	Material	Permeation	Thickness (mm)	Standard
Gloves	Nitrile rubber (NBR)	6 (> 480 minutes)	0,11	EN 374

#### Eye protection:

Safety glasses (to European standard EN 166 or equivalent)

#### Skin and body protection:

In case of possible repeated skin contact wear protective clothing

#### **Respiratory protection:**

Where excessive dust may result, wear approved mask. Type P1 (EN 143)

Safety Data Sheet

according to Regulation (EU) 2015/830

#### SECTION 9: Physical and chemical properties

#### 9.1. Information on basic physical and chemical properties

3.1. Information on basic physical and chemical properties		
Physical state	: Solid : Powder.	
Appearance		
Colour	: White to slightly yellow.	
Odour	: Characteristic. Weak.	
Odour threshold	: No data available	
рН	: No data available	
Relative evaporation rate (butylacetate=1)	: No data available	
Melting point	: No data available	
Freezing point	: No data available	
Boiling point	: No data available	
Flash point	: No data available	
Auto-ignition temperature	: No data available	
Decomposition temperature	: No data available	
Flammability (solid, gas)	: No data available	
Vapour pressure	: No data available	
Relative vapour density at 20 °C	: No data available	
Relative density	: No data available	
Solubility	: Readily soluble in water.	
Log Pow	: No data available	
Viscosity, kinematic	: No data available	
Viscosity, dynamic	: No data available	
Explosive properties	: No data available	
Oxidising properties	: No data available	
Explosive limits	: No data available	

#### 9.2. Other information

No additional information available

#### SECTION 10: Stability and reactivity

#### 10.1. Reactivity

Stable under normal conditions of storage, handling and use.

#### 10.2. Chemical stability

Stable under normal conditions.

#### 10.3. Possibility of hazardous reactions

No additional information available

#### **10.4.** Conditions to avoid

Moisture.

#### 10.5. Incompatible materials

Strong oxidizers.

#### 10.6. Hazardous decomposition products

Thermal decomposition generates :

- COx

Safety Data Sheet

according to Regulation (EU) 2015/830

- NOx

- SOx.

#### **SECTION 11: Toxicological information**

#### 11.1. Information on toxicological effects

Acute toxicity	: Not classified
Skin corrosion/irritation	: Not classified
Serious eye damage/irritation	: Not classified
Respiratory or skin sensitisation	: Not classified
Germ cell mutagenicity	: Not classified
Carcinogenicity	: Not classified
Reproductive toxicity	: Not classified
STOT-single exposure	: Not classified
STOT-repeated exposure	: Not classified
Aspiration hazard	: Not classified

#### **SECTION 12: Ecological information**

#### 12.1. Toxicity

No additional information available

#### 12.2. Persistence and degradability

No additional information available

#### 12.3. Bioaccumulative potential

No additional information available

#### 12.4. Mobility in soil

No additional information available

#### 12.5. Results of PBT and vPvB assessment

No additional information available

#### 12.6. Other adverse effects

Additional information

: Prevent entry to sewers and public waters. Avoid release to the environment

#### **SECTION 13: Disposal considerations**

#### 13.1. Waste treatment methods

Waste treatment methods

: Dispose in a safe manner in accordance with local/national regulations. Avoid release to the environment.

#### **SECTION 14: Transport information**

In accordance with ADR / IATA / IMDG

ADR	IMDG	ΙΑΤΑ
14.1. UN number		
Not regulated	Not regulated	Not regulated
17/10/2018 Duchefa Biochemie B.V.	EN (English)	SDS Ref.: L1704 5/7

Safety Data Sheet

according to Regulation (EU) 2015/830

ADR	IMDG	ΙΑΤΑ			
14.2. UN proper shipping name	14.2. UN proper shipping name				
Not regulated	Not regulated	Not regulated			
Not regulated	Not regulated	Not regulated			
14.3. Transport hazard class(es)					
Not regulated	Not regulated	Not regulated			
Not regulated	Not regulated	Not regulated			
14.4. Packing group					
Not regulated	Not regulated	Not regulated			
14.5. Environmental hazards					
Not regulated	Not regulated	Not regulated			
No supplementary information available					

#### 14.6. Special precautions for user

#### . Overland transport

Not regulated

#### - Transport by sea

Not regulated

#### - Air transport

Not regulated

#### 14.7. Transport in bulk according to Annex II of Marpol and the IBC Code

Not applicable

#### SECTION 15: Regulatory information

### 15.1. Safety, health and environmental regulations/legislation specific for the substance or mixture

#### 15.1.1. EU-Regulations

Contains no REACH substances with Annex XVII restrictions Contains no substance on the REACH candidate list Contains no REACH Annex XIV substances

#### 15.1.2. National regulations

Ensure all national/local regulations are observed.

#### Germany

Reference to AwSV

12th Ordinance Implementing the Federal Immission Control Act - 12.BImSchV

#### Netherlands

SZW-lijst van kankerverwekkende stoffen	: Yeast extract is listed
SZW-lijst van mutagene stoffen	: Yeast extract is listed
NIET-limitatieve lijst van voor de voortplanting giftige stoffen – Borstvoeding	: None of the components are listed
NIET-limitatieve lijst van voor de voortplanting giftige stoffen – Vruchtbaarheid	: None of the components are listed
NIET-limitatieve lijst van voor de voortplanting giftige stoffen – Ontwikkeling	: None of the components are listed

: Water hazard class (WGK) 3, severe hazard to water (Classification according to AwSV, Annex 1)

: Is not subject of the 12. BImSchV

(Hazardous Incident Ordinance)

Safety Data Sheet

according to Regulation (EU) 2015/830

#### 15.2. Chemical safety assessment

No additional information available

#### SECTION 16: Other information

#### Indication of changes:

	0		
8.2	Hand protection	Modified	Specified material, thickness, et
			cetera of gloves

Abbreviations and acronyms:

ATE	Acute Toxicity Estimate
ADR	European Agreement concerning the International Carriage of Dangerous Goods by Road
BCF	Bioconcentration factor
CLP	Classification Labelling Packaging Regulation; Regulation (EC) No 1272/2008
DPD	Dangerous Preparations Directive 1999/45/EC
DSD	Dangerous Substances Directive 67/548/EEC
IATA	International Air Transport Association
IMDG	International Maritime Dangerous Goods
LC50	Median lethal concentration
LD50	Median lethal dose
LOAEL	Lowest Observed Adverse Effect Level
NOAEC	No-Observed Adverse Effect Concentration
PBT	Persistent Bioaccumulative Toxic
REACH	Registration, Evaluation, Authorisation and Restriction of Chemicals Regulation (EC) No 1907/2006
SDS	Safety Data Sheet
Data sources	: REGULATION (EC) № 1272/2008 OF THE EUROPEAN PARLIAMENT

REGULATION (EC) No 1272/2008 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL of 16 December 2008 on classification, labelling and packaging of substances and mixtures, amending and repealing Directives 67/548/EEC and 1999/45/EC, and amending Regulation (EC) No 1907/2006. TNO (Netherlands Organisation for Applied Scientific Research).

#### SDS Biochemicals version 2018

This information is based on our current knowledge and is intended to describe the product for the purposes of health, safety and environmental requirements only. It should not therefore be construed as guaranteeing any specific property of the product