

L1717**Luria Broth Base, Miller**Medium: **General bacterial medium**Purpose: **Cultivation of bacteria****COMPOSITION OF MEDIA
L1717: LURIA BROTH BASE, MILLER**

COMPOUND	GRAM/LITER
Sodium chloride	0.5
Tryptone	10.0
Yeast Extract	5.0

METHOD

- Dissolve 16.5 grams of ingredients in distilled water and adjust volume to 1000 ml.
- Autoclave the solution (121 °C, 15 psi, 15 minutes).
- Allow medium to cool down.

L1717 LURIA BROTH BASE, MILLER

L 1717.0100	100 g
L 1717.0500	500 g
L 1717.1000	1 kg

For prepared and ready to use plates of this medium contact:
 Tritium Microbiologie Tel : 040-2051615
 Rooijakkersstraat 6 Fax : 040-2051395
 5652 BB Eindhoven Email : info@tritium-microbiologie.nl
 The Netherlands