

FICHE TECHNIQUE



M0132: MICONAZOLE NITRATE

CAS number 22832-87-7

Molecular weight C₁₈H₁₄Cl₄N₂O.HNO₃ 479.1

Storage condition

- Store at room temperature
- Soluble in propyleneglycol

Hazard statements

- **H302** Harmful if swallowed
- **H410** Very toxic to aquatic life with long lasting effects

Precautionary statements

- **P273** Avoid release to the environment
- **P280** Wear protective gloves/protective clothing/eye protection/face protection
- **P301** IF SWALLOWED:
- **P312** Call a POISON CENTER or doctor/physician if you feel unwell

M0132.0001/ 1 g

M0132.0005/ 5 g