

FICHE TECHNIQUE



M0135 : NEOMYCIN SULPHATE

CAS number 1405-10-3

Molecular weight $C_{23}H_{46}N_6O_{13} \cdot 3H_2SO_4 = 908.9$

Assay > 680 micro g / mg

Storage condition

- Store at room temperature
- Soluble in water (500 g/l / 20 °C)

Hazard statements

- **H315** Causes skin irritation
- **H317** May cause an allergic skin reaction
- **H319** Causes serious eye irritation
- **H334** May cause allergy or asthma symptoms or breathing difficulties if inhaled
- **H361** Suspected of damaging fertility or the unborn child

Precautionary statements

- **P261** Avoid breathing dust/fume/gas/mist/vapours/spray
- **P280** Wear protective gloves/protective clothing/eye protection/face protection
- **P305** IF IN EYES:
- **P311** Call a POISON CENTER or doctor/physician
- **P342** If experiencing respiratory symptoms:
- **P351** Rinse continuously with water for several minutes

M0135.0025/ 25 g

M0135.0100/ 100 g

