

FICHE TECHNIQUE



M0172 : MINOCYCLINE HCL

CAS number 13614-98-7

Molecular weight C₂₃H₂₇N₃O₇.HCl = 493.9

Assay > 96 %

Storage condition

- Store at 2-8 °C
- Slightly soluble in water
- Protect from light

Hazard statements

- **H315** Causes skin irritation
- **H319** Causes serious eye irritation
- **H335** May cause respiratory irritation
- **H360** May damage fertility or the unborn child
- **H362** May cause harm to breast-fed children

Precautionary statements

- **P201** Obtain special instructions before use
- **P263** Avoid contact during pregnancy/while nursing
- **P280** Wear protective gloves/protective clothing/eye protection/face protection
- **P308** IF exposed or concerned:
- **P313** Get medical advice/attention

M0172.0001/ 1 g