

FICHE TECHNIQUE



M0222 : MURASHIGE & SKOOG MEDIUM INCLUDING VITAMINS

Original concentration, (1962)

Micro Elements

mg/l

µM

CoCl₂.6H₂O

0.025

0.11

CuSO₄.5H₂O

0.025

0.10

FeNaEDTA

36.70

100.00

H₃BO₃

6.20

100.27

KI

0.83

5.00

MnSO₄.H₂O

16.90

100.00

Na₂MoO₄.2H₂O

0.25

1.03

ZnSO₄.7H₂O

8.60

29.91

Macro Elements

mg/l

mM

CaCl₂

332.02

2.99

KH₂PO₄

170.00

1.25
KNO₃
1900.00
18.79
MgSO₄
180.54
1.50
NH₄NO₃
1650.00
20.61

Vitamins

mg/l

μM

Glycine

2.00

26.64

myo-Inositol

100.00

554.94

Nicotinic acid

0.50

4.06

Pyridoxine HCl

0.50

2.43

Thiamine HCl

0.10

0.30

Total concentration Micro and Macro elements including vitamins: 4405.19 mg/l
Murashige T. and Skoog F., Physiol. Plant, 15, 473 (1962).

Hazard statements

- **H315** Causes skin irritation
- **H319** Causes serious eye irritation
- **H335** May cause respiratory irritation

Precautionary statements

- **P261** Avoid breathing dust/fume/gas/mist/vapours/spray
- **P302** IF ON SKIN:
- **P305** IF IN EYES:
- **P351** Rinse continuously with water for several minutes
- **P352** Wash with soap and water

M0222.0001/ 1 l
M0222.0005/ 5 l
M0222.0010/ 10 l
M0222.0025/ 25 l
M0222.0050/ 50 l
M0222.0075/ 75 l
M0222.0100/ 100 l

Parc Euromédecine -20, rue Robert Koch - BP 27253 34086 MONTPELLIER Cedex 4 - Tel : +33(0)4 67 10 59 00
RCS 461 800 559 B Montpellier - SIRET 461 800 559 00035 - CODE T.V.A CEE. : FR 62 461 800 559

contact@groupe-meridis.fr - www.groupe-meridis.fr