

FICHE TECHNIQUE



M0232 : MURASHIGE & SKOOG MEDIUM MOD.NO. 1A

Modification No.1 : 1/2 concentration Macro elem

Micro Elements

mg/l

μM

CoCl₂.6H₂O

0.025

0.11

CuSO₄.5H₂O

0.025

0.10

FeNaEDTA

36.70

100.00

H₃BO₃

6.20

100.27

KI

0.83

5.00

MnSO₄.H₂O

16.90

100.00

Na₂MoO₄.2H₂O

0.25

1.03

ZnSO₄.7H₂O

8.60

29.91

Macro Elements

mg/l

mM

CaCl₂

166.00

1.50

KH₂PO₄

85.00

0.63

KNO3

950.00

9.40

MgSO4

87.86

0.73

NH4NO3

825

10.3

Murashige T. and Skoog F., *Physiol. Plant*, 15, 473 (1962).

Hazard statements

- **H315** Causes skin irritation
- **H319** Causes serious eye irritation
- **H335** May cause respiratory irritation
- **H412** Harmful to aquatic life with long lasting effects

Precautionary statements

- **P261** Avoid breathing dust/fume/gas/mist/vapours/spray
- **P273** Avoid release to the environment
- **P302** IF ON SKIN:
- **P305** IF IN EYES:
- **P351** Rinse continuously with water for several minutes
- **P352** Wash with soap and water

M0232.0001/ 1 l

M0232.0010/ 10 l

M0232.0050/ 50 l

Parc Euromédecine -20, rue Robert Koch - BP 27253 34086 MONTPELLIER Cedex 4 - Tel : +33(0)4 67 10 59 00

RCS 461 800 559 B Montpellier - SIRET 461 800 559 00035 - CODE T.V.A CEE. : FR 62 461 800 559

contact@groupe-meridis.fr - www.groupe-meridis.fr