

Murashige & Skoog Medium mod. No.3B (Micro and Macro elements incl. vitamins, 1/2 conc. NH₄NO₃ and KNO₃)

Safety Data Sheet

according to Regulation (EU) 2015/830

Date of issue: 11/11/2016 Revision date 31/07/2018

M0237

Version: 2.0

SECTION 1: Identification of the substance/mixture and of the company/undertaking

1.1. Product identifier

Product form : Mixture
Trade name : Murashige & Skoog Medium mod. No.3B
(Micro and Macro elements incl. vitamins,
1/2 conc. NH₄NO₃ and KNO₃)
Product code : M0237
Product group : Raw material

1.2. Relevant identified uses of the substance or mixture and uses advised against

1.2.1. Relevant identified uses

Main use category : Professional use
Industrial/Professional use spec : For professional use only. Duchefa Biochemie B.V. products are intended only for "in vitro laboratory" research purposes.

1.2.2. Uses advised against

No additional information available

1.3. Details of the supplier of the safety data sheet

Manufacturer

Duchefa Biochemie B.V.
A. Hofmanweg 71
2031 BH Haarlem - The Netherlands
T +31(0)23-5319093 - F +31(0)23-5318027
info@duchefa.nl

1.4. Emergency telephone number

Emergency number : Supplier contact information:
+31(0)23-5319093 (M-F 09:00-17:00)
+31(0)6-30109355 (outside office hours)

Organisation/Company	Address	Comment
World Health Organization world directory of poison centres	http://apps.who.int/poisoncentres/	Consult website for a local poison centre

SECTION 2: Hazards identification

2.1. Classification of the substance or mixture

Classification according to Regulation (EC) No. 1272/2008 [CLP]

Skin corrosion/irritation, Category 2 H315
Serious eye damage/eye irritation, Category 2 H319
Specific target organ toxicity — Single exposure, Category 3, Respiratory tract irritation H335
Hazardous to the aquatic environment — Chronic Hazard, Category 3 H412

Full text of H statements : see section 16

Murashige & Skoog Medium mod. No.3B (Micro and Macro elements incl. vitamins, 1/2 conc. NH₄NO₃ and KNO₃)

M0237

Safety Data Sheet

according to Regulation (EU) 2015/830

Adverse physicochemical, human health and environmental effects

No additional information available

2.2. Label elements

Labelling according to Regulation (EC) No. 1272/2008 [CLP]

Hazard pictograms (CLP)



GHS07

Signal word (CLP)

: Warning

Hazard statements (CLP)

: H315 - Causes skin irritation.
H319 - Causes serious eye irritation.
H335 - May cause respiratory irritation.
H412 - Harmful to aquatic life with long lasting effects.

Precautionary statements (CLP)

: P261 - Avoid breathing dust.
P273 - Avoid release to the environment
P302+P352 - IF ON SKIN: Wash with plenty of water.
P305+P351 - IF IN EYES: Rinse cautiously with water for several minutes.

Extra phrases

: Based on research by TNO in Rijswijk (The Netherlands), commissioned by Duchefa Biochemie B.V. in Haarlem, the medium has no oxidising or explosive properties. As such the substance is not classified as oxidizing (H272, GHS03).

2.3. Other hazards

No additional information available

SECTION 3: Composition/information on ingredients

3.1. Substances

Not applicable

3.2. Mixtures

Name	Product identifier	%	Classification according to Regulation (EC) No. 1272/2008 [CLP]
Potassium nitrate	(CAS-No.) 7757-79-1 (EC-No.) 231-818-8	36,1191	Ox. Sol. 2, H272
Ammonium nitrate	(CAS-No.) 6484-52-2 (EC-No.) 229-347-8	31,3666	Ox. Sol. 3, H272
Calcium chloride	(CAS-No.) 10043-52-4 (EC-No.) 233-140-8 (EC Index-No.) 017-013-00-2	12,6234	Eye Irrit. 2, H319
Magnesium sulfate anhydrous	(CAS-No.) 7487-88-9 (EC-No.) 231-298-2	6,8641	Not classified
Potassium dihydrogenphosphate	(CAS-No.) 7778-77-0 (EC-No.) 231-913-4	6,4634	Not classified

Murashige & Skoog Medium mod. No.3B (Micro and Macro elements incl. vitamins, 1/2 conc. NH₄NO₃ and KNO₃)

M0237

Safety Data Sheet

according to Regulation (EU) 2015/830

Myo-Inositol	(CAS-No.) 87-89-8 (EC-No.) 201-781-2	3,802	Not classified
Ethylenediaminetetraacetate ferric sodium	(CAS-No.) 15708-41-5 (EC-No.) 239-802-2	1,3953	Not classified
Manganese sulfate monohydrate	(CAS-No.) 10034-96-5 (EC-No.) 232-089-9 (EC Index-No.) 025-003-00-4	0,6425	Eye Dam. 1, H318 STOT RE 2, H373 Aquatic Chronic 2, H411
Zinc sulphate	(CAS-No.) 7446-20-0 (EC-No.) 231-793-3 (EC Index-No.) 030-006-00-9	0,327	Acute Tox. 4 (Oral), H302 Eye Dam. 1, H318 Aquatic Acute 1, H400 Aquatic Chronic 1, H410
Boric acid substance listed as REACH Candidate	(CAS-No.) 10043-35-3 (EC-No.) 233-139-2 (EC Index-No.) 005-007-00-2	0,2357	Repr. 1B, H360FD
Glycine	(CAS-No.) 56-40-6 (EC-No.) 200-272-2	0,076	Not classified
Potassium iodide	(CAS-No.) 7681-11-0 (EC-No.) 231-659-4	0,0316	Acute Tox. 4 (Oral), H302 Skin Irrit. 2, H315 Eye Irrit. 2, H319
Pyridoxine HCl (Vitamine B6)	(CAS-No.) 58-56-0 (EC-No.) 200-386-2	0,019	Eye Dam. 1, H318
Nicotinic acid	(CAS-No.) 59-67-6 (EC-No.) 200-441-0	0,019	Eye Irrit. 2, H319
Sodium molybdate dihydrate	(CAS-No.) 10102-40-6 (EC-No.) 231-551-7	0,0095	Not classified
Thiamine hydrochloride	(CAS-No.) 67-03-8 (EC-No.) 200-641-8	0,0038	Not classified
Cobalt(II) chloride substance listed as REACH Candidate (Cobalt dichloride)	(CAS-No.) 7646-79-9 (EC-No.) 231-589-4 (EC Index-No.) 027-004-00-5	0,001	Acute Tox. 4 (Oral), H302 Resp. Sens. 1, H334 Skin Sens. 1, H317 Muta. 2, H341 Carc. 1B, H350i Repr. 1B, H360F Aquatic Acute 1, H400 Aquatic Chronic 1, H410
Copper sulphate	(CAS-No.) 7758-99-8 (EC-No.) 231-847-6 (EC Index-No.) 029-004-00-0	0,001	Acute Tox. 4 (Oral), H302 Skin Irrit. 2, H315 Eye Irrit. 2, H319 Aquatic Acute 1, H400 Aquatic Chronic 1, H410

Specific concentration limits:

Name	Product identifier	Specific concentration limits
Boric acid	(CAS-No.) 10043-35-3 (EC-No.) 233-139-2 (EC Index-No.) 005-007-00-2	(C >= 5,5) Repr. 1B, H360FD
Cobalt(II) chloride	(CAS-No.) 7646-79-9 (EC-No.) 231-589-4 (EC Index-No.) 027-004-00-5	(C >= 0,01) Carc. 1B, H350i

Full text of H-statements: see section 16

Murashige & Skoog Medium mod. No.3B (Micro and Macro elements incl. vitamins, 1/2 conc. NH₄NO₃ and KNO₃)

M0237

Safety Data Sheet

according to Regulation (EU) 2015/830

SECTION 4: First aid measures

4.1. Description of first aid measures

First-aid measures general	: Seek medical attention if ill effect develops.
First-aid measures after inhalation	: Remove victim to fresh air.
First-aid measures after skin contact	: Wash skin with mild soap and water.
First-aid measures after eye contact	: Rinse with water.
First-aid measures after ingestion	: Rinse mouth.

4.2. Most important symptoms and effects, both acute and delayed

Symptoms/effects after skin contact	: Causes skin irritation.
Symptoms/effects after eye contact	: Redness, pain.

4.3. Indication of any immediate medical attention and special treatment needed

No additional information available

SECTION 5: Firefighting measures

5.1. Extinguishing media

Suitable extinguishing media	: Alcohol resistant foam. Dry chemical powder. Carbon dioxide (CO ₂). Water spray.
------------------------------	---

5.2. Special hazards arising from the substance or mixture

Hazardous decomposition products in case of fire	: Under fire conditions, hazardous fumes will be present: - CO _x - NO _x - SO _x .
--	--

5.3. Advice for firefighters

Firefighting instructions	: Prevent fire fighting water from entering the environment.
Protection during firefighting	: Wear proper protective equipment.

SECTION 6: Accidental release measures

6.1. Personal precautions, protective equipment and emergency procedures

General measures	: Avoid raising powdered materials into airborne dust.
------------------	--

6.1.1. For non-emergency personnel

Emergency procedures	: Wear suitable protective clothing.
----------------------	--------------------------------------

6.1.2. For emergency responders

No additional information available

6.2. Environmental precautions

Prevent entry to sewers and public waters.

Murashige & Skoog Medium mod. No.3B (Micro and Macro elements incl. vitamins, 1/2 conc. NH₄NO₃ and KNO₃)

M0237

Safety Data Sheet

according to Regulation (EU) 2015/830

6.3. Methods and material for containment and cleaning up

Methods for cleaning up : Sweep up dry powder and dispose properly.

6.4. Reference to other sections

No additional information available

SECTION 7: Handling and storage

7.1. Precautions for safe handling

Precautions for safe handling : Avoid dust formation. Handle in accordance with good industrial hygiene and safety procedures.

7.2. Conditions for safe storage, including any incompatibilities

Storage conditions : Store at room temperature
Store in dry, well-ventilated area
Hygroscopic
Keep container tightly closed and dry.

7.3. Specific end use(s)

For professional use only. Duchefa Biochemie B.V. products are intended only for "in vitro laboratory" research purposes.

SECTION 8: Exposure controls/personal protection

8.1. Control parameters

No additional information available

8.2. Exposure controls

Hand protection:

Type	Material	Permeation	Thickness (mm)	Standard
Gloves	Nitrile rubber (NBR)	6 (> 480 minutes)	0,11	EN 374

Eye protection:

Safety glasses (to European standard EN 166 or equivalent)

Skin and body protection:

In case of possible repeated skin contact wear protective clothing

Respiratory protection:

Where excessive dust may result, wear approved mask. Type P2 (EN 143)

SECTION 9: Physical and chemical properties

9.1. Information on basic physical and chemical properties

Physical state : Solid
Appearance : Powder.
Colour : White to slightly yellow.
Odour : Characteristic. Weak.
Odour threshold : No data available

Murashige & Skoog Medium mod. No.3B (Micro and Macro elements incl. vitamins, 1/2 conc. NH₄NO₃ and KNO₃)

M0237

Safety Data Sheet

according to Regulation (EU) 2015/830

pH	: No data available
Relative evaporation rate (butylacetate=1)	: No data available
Melting point	: No data available
Freezing point	: No data available
Boiling point	: No data available
Flash point	: No data available
Auto-ignition temperature	: No data available
Decomposition temperature	: No data available
Flammability (solid, gas)	: No data available
Vapour pressure	: No data available
Relative vapour density at 20 °C	: No data available
Relative density	: No data available
Solubility	: Readily soluble in water.
Log Pow	: No data available
Viscosity, kinematic	: No data available
Viscosity, dynamic	: No data available
Explosive properties	: No data available
Oxidising properties	: No data available
Explosive limits	: No data available

9.2. Other information

No additional information available

SECTION 10: Stability and reactivity

10.1. Reactivity

Stable under normal conditions of storage, handling and use.

10.2. Chemical stability

Stable under normal conditions.

10.3. Possibility of hazardous reactions

No additional information available

10.4. Conditions to avoid

Moisture.

10.5. Incompatible materials

Strong oxidizers.

10.6. Hazardous decomposition products

Thermal decomposition generates :

- CO_x
- NO_x
- SO_x
- Phosphorus oxides.

Murashige & Skoog Medium mod. No.3B (Micro and Macro elements incl. vitamins, 1/2 conc. NH₄NO₃ and KNO₃)

M0237

Safety Data Sheet

according to Regulation (EU) 2015/830

SECTION 11: Toxicological information

11.1. Information on toxicological effects

Acute toxicity	: Not classified
Skin corrosion/irritation	: Causes skin irritation.
Serious eye damage/irritation	: Causes serious eye irritation.
Respiratory or skin sensitisation	: Not classified
Germ cell mutagenicity	: Not classified
Carcinogenicity	: Not classified
Reproductive toxicity	: Not classified
STOT-single exposure	: May cause respiratory irritation.
STOT-repeated exposure	: Not classified
Aspiration hazard	: Not classified

SECTION 12: Ecological information

12.1. Toxicity

No additional information available

12.2. Persistence and degradability

No additional information available

12.3. Bioaccumulative potential

No additional information available

12.4. Mobility in soil

No additional information available

12.5. Results of PBT and vPvB assessment

Murashige & Skoog Medium mod. No.3B (Micro and Macro elements incl. vitamins, 1/2 conc. NH₄NO₃ and KNO₃)	
Component	
Boric acid (10043-35-3)	This substance/mixture does not meet the PBT criteria of REACH regulation, annex XIII This substance/mixture does not meet the vPvB criteria of REACH regulation, annex XIII

12.6. Other adverse effects

Additional information : Prevent entry to sewers and public waters. Avoid release to the environment

SECTION 13: Disposal considerations

13.1. Waste treatment methods

Waste treatment methods : Dispose in a safe manner in accordance with local/national regulations.
Avoid release to the environment.

Murashige & Skoog Medium mod. No.3B (Micro and Macro elements incl. vitamins, 1/2 conc. NH₄NO₃ and KNO₃)

M0237

Safety Data Sheet

according to Regulation (EU) 2015/830

SECTION 14: Transport information

In accordance with ADR / IATA / IMDG

ADR	IMDG	IATA
14.1. UN number		
Not regulated	Not regulated	Not regulated
14.2. UN proper shipping name		
Not regulated	Not regulated	Not regulated
Not regulated	Not regulated	Not regulated
14.3. Transport hazard class(es)		
Not regulated	Not regulated	Not regulated
Not regulated	Not regulated	Not regulated
14.4. Packing group		
Not regulated	Not regulated	Not regulated
14.5. Environmental hazards		
Not regulated	Not regulated	Not regulated
No supplementary information available		

14.6. Special precautions for user

Overland transport

Not regulated

- Transport by sea

Not regulated

- Air transport

Not regulated

14.7. Transport in bulk according to Annex II of Marpol and the IBC Code

Not applicable

SECTION 15: Regulatory information

15.1. Safety, health and environmental regulations/legislation specific for the substance or mixture

15.1.1. EU-Regulations

Contains no REACH substances with Annex XVII restrictions

Contains a substance on the REACH candidate list in concentration $\geq 0.1\%$ or with a lower specific limit: Boric acid (EC 233-139-2, CAS 10043-35-3)

Contains no REACH Annex XIV substances

15.1.2. National regulations

Ensure all national/local regulations are observed.

Germany

Reference to AwSV

: Water hazard class (WGK) 2,
significant hazard to water
(Classification according to AwSV,
Annex 1)

Murashige & Skoog Medium mod. No.3B (Micro and Macro elements incl. vitamins, 1/2 conc. NH₄NO₃ and KNO₃)

M0237

Safety Data Sheet

according to Regulation (EU) 2015/830

12th Ordinance Implementing the Federal Immission Control Act -
12.BImSchV

: Is not subject of the 12. BImSchV
(Hazardous Incident Ordinance)

Netherlands

SZW-lijst van kankerverwekkende stoffen

: Manganese sulfate
monohydrate, Cobalt(II) chloride are
listed

SZW-lijst van mutagene stoffen

: Manganese sulfate monohydrate is
listed

NIET-limitatieve lijst van voor de voortplanting giftige stoffen – Borstvoeding

: None of the components are listed

NIET-limitatieve lijst van voor de voortplanting giftige stoffen –
Vruchtbaarheid

: Boric acid, Cobalt(II) chloride are
listed

NIET-limitatieve lijst van voor de voortplanting giftige stoffen – Ontwikkeling

: Boric acid is listed

Denmark

Recommendations Danish Regulation

: Young people below the age of 18
years are not allowed to use the
product

Pregnant/breastfeeding women
working with the product must not be
in direct contact with the product

The requirements from the Danish
Working Environment Authorities
regarding work with carcinogens must
be followed during use and disposal

15.2. Chemical safety assessment

No additional information available

SECTION 16: Other information

Indication of changes:

2.2	Signal word (CLP)	Added	Warning
2.2	Hazard pictograms (CLP)	Added	GHS07
2.2	Hazard statements (CLP)	Added	H315;H319;H335;H412
2.2	Precautionary statements (CLP)	Added	P261;P273;P302+P352;P305+P351+P338

Abbreviations and acronyms:

ATE	Acute Toxicity Estimate
ADR	European Agreement concerning the International Carriage of Dangerous Goods by Road
BCF	Bioconcentration factor
CLP	Classification Labelling Packaging Regulation; Regulation (EC) No 1272/2008
DPD	Dangerous Preparations Directive 1999/45/EC
DSD	Dangerous Substances Directive 67/548/EEC
IATA	International Air Transport Association
IMDG	International Maritime Dangerous Goods
LC50	Median lethal concentration
LD50	Median lethal dose
LOAEL	Lowest Observed Adverse Effect Level
NOAEC	No-Observed Adverse Effect Concentration
PBT	Persistent Bioaccumulative Toxic

Murashige & Skoog Medium mod. No.3B (Micro and Macro elements incl. vitamins, 1/2 conc. NH₄NO₃ and KNO₃)

M0237

Safety Data Sheet

according to Regulation (EU) 2015/830

REACH	Registration, Evaluation, Authorisation and Restriction of Chemicals Regulation (EC) No 1907/2006
SDS	Safety Data Sheet

Data sources : REGULATION (EC) No 1272/2008 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL of 16 December 2008 on classification, labelling and packaging of substances and mixtures, amending and repealing Directives 67/548/EEC and 1999/45/EC, and amending Regulation (EC) No 1907/2006. TNO (Netherlands Organisation for Applied Scientific Research).

Full text of H- and EUH-statements:

Acute Tox. 4 (Oral)	Acute toxicity (oral), Category 4
Aquatic Acute 1	Hazardous to the aquatic environment — Acute Hazard, Category 1
Aquatic Chronic 1	Hazardous to the aquatic environment — Chronic Hazard, Category 1
Aquatic Chronic 2	Hazardous to the aquatic environment — Chronic Hazard, Category 2
Aquatic Chronic 3	Hazardous to the aquatic environment — Chronic Hazard, Category 3
Carc. 1B	Carcinogenicity (inhalation) Category 1B
Eye Dam. 1	Serious eye damage/eye irritation, Category 1
Eye Irrit. 2	Serious eye damage/eye irritation, Category 2
Muta. 2	Germ cell mutagenicity, Category 2
Ox. Sol. 2	Oxidising Solids, Category 2
Ox. Sol. 3	Oxidising Solids, Category 3
Repr. 1B	Reproductive toxicity, Category 1B
Repr. 1B	Reproductive toxicity, Category 1B
Resp. Sens. 1	Respiratory sensitisation, Category 1
Skin Irrit. 2	Skin corrosion/irritation, Category 2
Skin Sens. 1	Skin sensitisation, Category 1
STOT RE 2	Specific target organ toxicity — Repeated exposure, Category 2
STOT SE 3	Specific target organ toxicity — Single exposure, Category 3, Respiratory tract irritation
H272	May intensify fire; oxidiser.
H302	Harmful if swallowed.
H315	Causes skin irritation.
H317	May cause an allergic skin reaction.
H318	Causes serious eye damage.
H319	Causes serious eye irritation.
H334	May cause allergy or asthma symptoms or breathing difficulties if inhaled.
H335	May cause respiratory irritation.
H341	Suspected of causing genetic defects.
H350i	May cause cancer by inhalation.
H360F	May damage fertility.
H360FD	May damage fertility. May damage the unborn child.
H373	May cause damage to organs through prolonged or repeated exposure.
H400	Very toxic to aquatic life.
H410	Very toxic to aquatic life with long lasting effects
H411	Toxic to aquatic life with long lasting effects.
H412	Harmful to aquatic life with long lasting effects.

SDS Biochemicals version 2018

**Murashige & Skoog Medium mod. No.3B
(Micro and Macro elements incl. vitamins,
1/2 conc. NH₄NO₃ and KNO₃)**

M0237

Safety Data Sheet

according to Regulation (EU) 2015/830

This information is based on our current knowledge and is intended to describe the product for the purposes of health, safety and environmental requirements only. It should not therefore be construed as guaranteeing any specific property of the product