

FICHE TECHNIQUE



M0301: MICRO SALT MIXTURE MS

Micro Elements

mg/l

μM

CoCl₂.6H₂O

0.025

0.11

CuSO₄.5H₂O

0.025

0.10

FeNaEDTA

36.70

100.00

H₃BO₃

6.20

100.27

KI

0.83

5.00

MnSO₄.H₂O

16.90

100.00

Na₂MoO₄.2H₂O

0.25

1.03

ZnSO₄.7H₂O

8.60

29.91

Macro Elements

mg/l

mM

KNO₃

930.47

9.20

Murashige T. and Skoog F., *Physiol. Plant*, 15, 473 (1962).

Hazard statements

- **H272** May intensify fire; oxidizer
- **H412** Harmful to aquatic life with long lasting effects

Precautionary statements

- **P210** Keep away from heat/sparks/open flames/hot surfaces
- **P273** Avoid release to the environment

M0301.0025/ 25 l

M0301.0050/ 50 l