

FICHE TECHNIQUE



M0303: MICRO SALT MIXTURE NITSCH

Micro Elements

mg/l

μ M

CuSO₄.5H₂O

0.025

FeNaEDTA

36.70

H₃BO₃

10.00

MnSO₄.H₂O

18.94

Na₂MoO₄.2H₂O

0.25

ZnSO₄.7H₂O

10.00

Macro Elements

mg/l

mM

KNO₃

424.85

Nitsch J.P. and Nitsch C., Haploid plants from pollen grains, Science 169, 85 (1969).

Nitsch J.P., Experimental androgenesis in Nicotiana, Phytomorphology 19, 389 (1969).

Hazard statements

- **H272** May intensify fire; oxidizer
- **H319** Causes serious eye irritation
- **H412** Harmful to aquatic life with long lasting effects

Precautionary statements

- **P210** Keep away from heat/sparks/open flames/hot surfaces
- **P273** Avoid release to the environment
- **P280** Wear protective gloves/protective clothing/eye protection/face protection
- **P305** IF IN EYES:
- **P351** Rinse continuously with water for several minutes

M0303.0025/

25 l