

## FICHE TECHNIQUE



### **M0304: MACRO SALT MIXTURE B5**

#### Macro Elements

mg/l

mM

CaCl<sub>2</sub>

113.23

1.02

KNO<sub>3</sub>

500.77

4.95

MgSO<sub>4</sub>

121.56

1.01

NaH<sub>2</sub>PO<sub>4</sub>

130.44

1.09

(NH<sub>4</sub>)<sub>2</sub>SO<sub>4</sub>

134.00

1.01

Gamborg O.L., Miller R.A., Ojima K., Nutrient requirement of suspensions cultures of soybean root cells. Exp. Cell Res., 50, 151 (1968).

#### Hazard statements

- **H272** May intensify fire; oxidizer
- **H319** Causes serious eye irritation
- **H373** May cause damage to organs through prolonged or repeated exposure
- **H412** Harmful to aquatic life with long lasting effects

#### Precautionary statements

- **P210** Keep away from heat/sparks/open flames/hot surfaces
- **P260** Do not breathe dust/fume/gas/mist/vapours/spray
- **P273** Avoid release to the environment
- **P280** Wear protective gloves/protective clothing/eye protection/face protection
- **P305** IF IN EYES:
- **P351** Rinse continuously with water for several minutes

M0304.0025/ 25 l