

FICHE TECHNIQUE



M0305: MACRO SALT MIXTURE MS

Macro Elements

mg/l

mM

CaCl₂

332.02

2.99

KH₂PO₄

170.00

1.25

KNO₃

969.53

9.59

MgSO₄

180.54

1.50

Murashige T. and Skoog F., *Physiol. Plant*, 15, 473 (1962).

Hazard statements

- **H272** May intensify fire; oxidizer
- **H319** Causes serious eye irritation

Precautionary statements

- **P210** Keep away from heat/sparks/open flames/hot surfaces
- **P280** Wear protective gloves/protective clothing/eye protection/face protection
- **P305** IF IN EYES:
- **P351** Rinse continuously with water for several minutes

M0305.0025/ 25 l

M0305.0050/ 50 l