

FICHE TECHNIQUE



M0306: MACRO SALT MIXTURE NITSCH

Macro Elements

mg/l

mM

CaCl₂

166.00

1.5

KH₂PO₄

68.00

0.5

KNO₃

525.92

5.2

MgSO₄

90.27

0.75

Nitsch J.P. and Nitsch C., Haploid plants from pollen grains, Science 169, 85 (1969).

Nitsch J.P., Experimental androgenesis in Nicotiana, Phytomorphology 19, 389 (1969).

Hazard statements

- **H272** May intensify fire; oxidizer
- **H319** Causes serious eye irritation

Precautionary statements

- **P210** Keep away from heat/sparks/open flames/hot surfaces
- **P280** Wear protective gloves/protective clothing/eye protection/face protection
- **P305** IF IN EYES:
- **P351** Rinse continuously with water for several minutes

M0306.0025/ 25 l