

FICHE TECHNIQUE

M0408: MCCOWN WOODY PLANT VITAMIN MIXTURE

Vitamins

mg/l

µM

Glycine

2.00

26.64

myo-Inositol

100.00

554.94

Nicotinic acid

0.50

4.06

Pyridoxine HCl

0.50

2.43

Thiamine HCl

1.00

2.96

Lloyd G. and McCown. Commercially-feasible micropropagation of mountain laurel, *Kalmia latifolia*, by use of shoot-tip culture. B., Int. Plant Prop. Soc. Proc. 30, 421 (1980).

M0408.0100/ >100 ml

M0408.0250/ >250 ml