

FICHE TECHNIQUE

M0409: MURASHIGE & SKOOG MEDIUM VITAMIN MIXTURE

Vitamins

mg/l

µM

Glycine

2.00

26.64

myo-Inositol

100.00

554.94

Nicotinic acid

0.50

4.06

Pyridoxine HCl

0.50

2.43

Thiamine HCl

0.10

0.30

Murashige T. and Skoog F., *Physiol. Plant*, 15, 473 (1962).

M0409.0100/ >100 ml

M0409.0250/ >250 ml