

FICHE TECHNIQUE



M0514: MANGANESE SULPHATE MONOHYDRATE

CAS number 10034-96-5

Molecular weight $\text{MnSO}_4 \cdot \text{H}_2\text{O}$ = 169.0

Assay > 98 %

Storage condition

- Store at room temperature
- Hygroscopic
- Soluble in water (750 g/l / 20 °C)

Hazard statements

- **H318** Causes serious eye damage
- **H373** May cause damage to organs through prolonged or repeated exposure
- **H411** Toxic to aquatic life with long lasting effects

Precautionary statements

- **P273** Avoid release to the environment
- **P280** Wear protective gloves/protective clothing/eye protection/face protection
- **P305** IF IN EYES:
- **P314** Get Medical advice/attention if you feel unwell
- **P351** Rinse continuously with water for several minutes

M0514.0250/ 250 g

M0514.0500/ 500 g

M0514.1000/ 1 kg