

Magnesium chloride hexahydrate

Safety Data Sheet

according to Regulation (EU) 2015/830

Date of issue: 14/09/2011 Revision date 15/06/2018

M0533

Version: 2.0

SECTION 1: Identification of the substance/mixture and of the company/undertaking

1.1. Product identifier

Product form : Substance
 Trade name : Magnesium chloride hexahydrate
 EC-No. : 232-094-6
 CAS-No. : 7791-18-6
 Product code : M0533
 Formula : MgCl₂ · 6H₂O

1.2. Relevant identified uses of the substance or mixture and uses advised against

1.2.1. Relevant identified uses

Main use category : Professional use
 Industrial/Professional use spec : For professional use only. Duchefa Biochemie B.V. products are intended only for "in vitro laboratory" research purposes.

1.2.2. Uses advised against

No additional information available

1.3. Details of the supplier of the safety data sheet

Supplier

Duchefa Biochemie B.V.
 A. Hofmanweg 71
 2031 BH Haarlem - The Netherlands
 T +31(0)23-5319093 - F +31(0)23-5318027
info@duchefa.nl

1.4. Emergency telephone number

Emergency number : Supplier contact information:
 +31(0)23-5319093 (M-F 09:00-17:00)
 +31(0)6-30109355 (outside office hours)

Organisation/Company	Address	Comment
World Health Organization world directory of poison centres	http://apps.who.int/poisoncentres/	Consult website for a local poison centre

SECTION 2: Hazards identification

2.1. Classification of the substance or mixture

Classification according to Regulation (EC) No. 1272/2008 [CLP]

Not classified

Adverse physicochemical, human health and environmental effects

No additional information available

2.2. Label elements

Labelling according to Regulation (EC) No. 1272/2008 [CLP]

No labelling applicable

2.3. Other hazards

No additional information available

SECTION 3: Composition/information on ingredients

3.1. Substances

Name : Magnesium chloride hexahydrate

CAS-No. : 7791-18-6

EC-No. : 232-094-6

Name	Product identifier	%
Magnesium chloride hexahydrate	(CAS-No.) 7791-18-6 (EC-No.) 232-094-6	> 98

Full text of H-statements: see section 16

3.2. Mixtures

Not applicable

SECTION 4: First aid measures

4.1. Description of first aid measures

First-aid measures after inhalation : Assure fresh air breathing.

First-aid measures after skin contact : Remove affected clothing and wash all exposed skin area with mild soap and water, followed by warm water rinse
Seek medical attention if ill effect or irritation develops.

First-aid measures after eye contact : Rinse with water
Consult an ophtalmologist if irritation persists.

First-aid measures after ingestion : Rinse mouth
Give water to drink.

4.2. Most important symptoms and effects, both acute and delayed

Symptoms/effects : Diarrhea
Vomiting
Affects the nervous system.

4.3. Indication of any immediate medical attention and special treatment needed

No additional information available

SECTION 5: Firefighting measures

5.1. Extinguishing media

Suitable extinguishing media : Dry chemical, CO₂, or water spray or regular foam.

5.2. Special hazards arising from the substance or mixture

Hazardous decomposition products in case of fire : Under fire conditions, hazardous fumes will be present.

5.3. Advice for firefighters

Firefighting instructions : Prevent fire fighting water from entering the environment.

Protection during firefighting : Do not enter fire area without proper protective equipment, including respiratory protection.

SECTION 6: Accidental release measures

6.1. Personal precautions, protective equipment and emergency procedures

General measures : Spill should be handled by trained cleaning personnel properly equipped with respiratory and eye protection.

6.1.1. For non-emergency personnel

Protective equipment : Wear recommended personal protective equipment.

Measures in case of dust release : Avoid dust formation.

6.1.2. For emergency responders

No additional information available

6.2. Environmental precautions

Prevent product from entering drains. Avoid release to the environment.

6.3. Methods and material for containment and cleaning up

Methods for cleaning up : Sweep up dry powder and dispose properly.

6.4. Reference to other sections

No additional information available

SECTION 7: Handling and storage

7.1. Precautions for safe handling

Precautions for safe handling : Ensure good ventilation of the work station. Handle in accordance with good industrial hygiene and safety procedures.

7.2. Conditions for safe storage, including any incompatibilities

Storage conditions : Store in a well-ventilated place. Keep container tightly closed.
Store at room temperature.

7.3. Specific end use(s)

For professional use only. Duchefa Biochemie B.V. products are intended only for "in vitro laboratory" research purposes.

SECTION 8: Exposure controls/personal protection

8.1. Control parameters

No additional information available

8.2. Exposure controls

Hand protection:

Type	Material	Permeation	Thickness (mm)	Standard
Gloves	Nitrile rubber (NBR)	6 (> 480 minutes)	0,11	EN 374

Eye protection:

Safety glasses (to European standard EN 166 or equivalent)

Skin and body protection:

Wear suitable protective clothing

Respiratory protection:

Dust production: dust mask with filter type P1

SECTION 9: Physical and chemical properties

9.1. Information on basic physical and chemical properties

Physical state	: Solid
Appearance	: Crystalline powder.
Molecular mass	: 203,3 g/mol
Colour	: White.
Odour	: Weak.
Odour threshold	: No data available
pH	: No data available
pH solution	: 4,5 - 7,5 (50 g/l)
Relative evaporation rate (butylacetate=1)	: No data available
Melting point	: 117 °C
Freezing point	: No data available
Boiling point	: No data available
Flash point	: No data available
Auto-ignition temperature	: No data available
Decomposition temperature	: No data available
Flammability (solid, gas)	: No data available
Vapour pressure	: No data available
Relative vapour density at 20 °C	: No data available
Relative density	: ≈ 1,57 (g/cm ³)
Solubility	: Water: 1670 g/l
Log Pow	: No data available
Viscosity, kinematic	: No data available
Viscosity, dynamic	: No data available
Explosive properties	: No data available
Oxidising properties	: No data available
Explosive limits	: No data available

9.2. Other information

No additional information available

SECTION 10: Stability and reactivity

10.1. Reactivity

Stable under normal conditions of storage, handling and use.

10.2. Chemical stability

Stable under normal conditions.

10.3. Possibility of hazardous reactions

No additional information available

10.4. Conditions to avoid

Moisture

Heat.

10.5. Incompatible materials

Strong oxidizers.

10.6. Hazardous decomposition products

When heated to decomposition, emits dangerous fumes.

SECTION 11: Toxicological information

11.1. Information on toxicological effects

Acute toxicity : Not classified

Magnesium chloride hexahydrate (7791-18-6)	
LD50 oral rat	8100 mg/kg
LD50 oral	7600 mg/kg (mouse)

Skin corrosion/irritation : Not classified

Serious eye damage/irritation : Not classified

Respiratory or skin sensitisation : Not classified

Germ cell mutagenicity : Not classified

Carcinogenicity : Not classified

Reproductive toxicity : Not classified

STOT-single exposure : Not classified

STOT-repeated exposure : Not classified

Aspiration hazard : Not classified

Other information : See actual entry in RTECS for complete information: OM2975000.

SECTION 12: Ecological information

12.1. Toxicity

Magnesium chloride hexahydrate (7791-18-6)	
LC50 fish 1	2120 mg/l Pimephales promelas
EC50 Daphnia 1	1400 mg/l
EC50 72h algae (1)	2200 mg/l

12.2. Persistence and degradability

No additional information available

12.3. Bioaccumulative potential

No additional information available

12.4. Mobility in soil

No additional information available

12.5. Results of PBT and vPvB assessment

No additional information available

12.6. Other adverse effects

Additional information : Avoid release to the environment. Prevent entry to sewers and public waters

SECTION 13: Disposal considerations**13.1. Waste treatment methods**

Waste treatment methods : Dispose in a safe manner in accordance with local/national regulations.

SECTION 14: Transport information

In accordance with ADR / IATA / IMDG

ADR	IMDG	IATA
14.1. UN number		
Not regulated	Not regulated	Not regulated
14.2. UN proper shipping name		
Not regulated	Not regulated	Not regulated
Not regulated	Not regulated	Not regulated
14.3. Transport hazard class(es)		
Not regulated	Not regulated	Not regulated
Not regulated	Not regulated	Not regulated
14.4. Packing group		
Not regulated	Not regulated	Not regulated
14.5. Environmental hazards		
Not regulated	Not regulated	Not regulated
No supplementary information available		

14.6. Special precautions for user**. Overland transport**

Not regulated

- Transport by sea

Not regulated

- Air transport

Not regulated

14.7. Transport in bulk according to Annex II of Marpol and the IBC Code

Not applicable

SECTION 15: Regulatory information**15.1. Safety, health and environmental regulations/legislation specific for the substance or mixture****15.1.1. EU-Regulations**

No REACH Annex XVII restrictions

Magnesium chloride hexahydrate is not on the REACH Candidate List

Magnesium chloride hexahydrate is not on the REACH Annex XIV List

15.1.2. National regulations

Ensure all national/local regulations are observed.

Germany

Magnesium chloride hexahydrate

M0533

Safety Data Sheet

according to Regulation (EU) 2015/830

Reference to AwSV	: Water hazard class (WGK) 1, low hazard to waters (Classification according to VwVwS, Annex 1 or 2; ID No. 259)
12th Ordinance Implementing the Federal Immission Control Act - 12.BImSchV	: Is not subject of the 12. BImSchV (Hazardous Incident Ordinance)
Netherlands	
SZW-lijst van kankerverwekkende stoffen	: Magnesium chloride hexahydrate is listed
SZW-lijst van mutagene stoffen	: Magnesium chloride hexahydrate is listed
NIET-limitatieve lijst van voor de voortplanting giftige stoffen – Borstvoeding	: The substance is not listed
NIET-limitatieve lijst van voor de voortplanting giftige stoffen – Vruchtbaarheid	: The substance is not listed
NIET-limitatieve lijst van voor de voortplanting giftige stoffen – Ontwikkeling	: The substance is not listed

15.2. Chemical safety assessment

No additional information available

SECTION 16: Other information

Indication of changes:

1.4	Emergency number	Modified	
8.2	Hand protection	Modified	Specified material, thickness, et cetera of gloves
11.1	Other information	Added	RTECS no
11.1	LD50 oral	Added	

Abbreviations and acronyms:

ATE	Acute Toxicity Estimate
ADR	European Agreement concerning the International Carriage of Dangerous Goods by Road
BCF	Bioconcentration factor
CLP	Classification Labelling Packaging Regulation; Regulation (EC) No 1272/2008
DPD	Dangerous Preparations Directive 1999/45/EC
DSD	Dangerous Substances Directive 67/548/EEC
IATA	International Air Transport Association
IMDG	International Maritime Dangerous Goods
LC50	Median lethal concentration
LD50	Median lethal dose
LOAEL	Lowest Observed Adverse Effect Level
NOAEC	No-Observed Adverse Effect Concentration
PBT	Persistent Bioaccumulative Toxic
REACH	Registration, Evaluation, Authorisation and Restriction of Chemicals Regulation (EC) No 1907/2006
SDS	Safety Data Sheet

Data sources : ECHA (European Chemicals Agency). Supplier. REGULATION (EC) No 1272/2008 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL of 16 December 2008 on classification, labelling and packaging of substances and mixtures, amending and repealing Directives 67/548/EEC and 1999/45/EC, and amending Regulation (EC) No 1907/2006.

SDS Biochemicals version 2018

Magnesium chloride hexahydrate

M0533

Safety Data Sheet

according to Regulation (EU) 2015/830

This information is based on our current knowledge and is intended to describe the product for the purposes of health, safety and environmental requirements only. It should not therefore be construed as guaranteeing any specific property of the product