

FICHE TECHNIQUE



M0921 : MALEIC HYDRAZIDE

CAS number 123-33-1

Molecular weight $C_4H_4N_2O_2 = 112.1$

Assay > 98 %

Storage condition

- Store dry at room temperature
- Soluble in 1 N NaOH, Insoluble in water
- 3,6-pyridazinediol

Hazard statements

- **H315** Causes skin irritation
- **H319** Causes serious eye irritation
- **H335** May cause respiratory irritation
- **H341** Suspected of causing genetic defects

Precautionary statements

- **P280** Wear protective gloves/protective clothing/eye protection/face protection
- **P302** IF ON SKIN:
- **P305** IF IN EYES:
- **P351** Rinse continuously with water for several minutes
- **P352** Wash with soap and water

M0921.0100/ 100 g

Parc Euromédecine -20, rue Robert Koch - BP 27253 34086 MONTPELLIER Cedex 4 - Tel : +33(0)4 67 10 59 00
RCS 461 800 559 B Montpellier - SIRET 461 800 559 00035 - CODE T.V.A CEE. : FR 62 461 800 559
contact@groupe-meridis.fr - www.groupe-meridis.fr