

Safety Data Sheet

M1315

Version: 2.0

according to Regulation (EU) 2015/830

Date of issue: 04/06/2010 Revision date: 21/11/2017 Supersedes 05/08/2011

SECTION 1: Identification of the substance/mixture and of the company/undertaking

1.1. **Product identifier**

Product form : Substance Trade name : (DL)-Malic acid

IUPAC name

EC-No. : 230-022-8 CAS-No. : 6915-15-7 Product code : M1315 Formula : C4H6O5

1.2. Relevant identified uses of the substance or mixture and uses advised against

1.2.1. Relevant identified uses

Main use category : Professional use

Industrial/Professional use spec : For professional use only. Duchefa Biochemie B.V. products are intended

only for "in vitro laboratory" research purposes.

1.2.2. Uses advised against

No additional information available

1.3. Details of the supplier of the safety data sheet

Duchefa Biochemie B.V. A. Hofmanweg 71 2031 BH Haarlem - The Netherlands T +31(0)23-5319093 - F +31(0)23-5318027 info@duchefa.nl

1.4. **Emergency telephone number**

Emergency number : Supplier contact information:

> +31(0)23-5319093 (M-F 09:00-17:00) +31(0)6-30109355 (outside office hours)

> > SDS Ref.: M1315

Organisation/Company	Address	Comment
World Health Organization world directory of poison centres	http://apps.who.int/poisoncentres/	Consult website for a local poison centre

SECTION 2: Hazards identification

2.1. Classification of the substance or mixture

Classification according to Regulation (EC) No. 1272/2008 [CLP]

Serious eye damage/eye irritation, Category 2 H319

Full text of hazard classes and H-statements : see section 16

Adverse physicochemical, human health and environmental effects

No additional information available

Safety Data Sheet

2.2. Label elements

Labelling according to Regulation (EC) No. 1272/2008 [CLP]

Hazard pictograms (CLP)



GHS07

Signal word (CLP) : Warning

Hazard statements (CLP) : H319 - Causes serious eye irritation.

Precautionary statements (CLP) : P280 - Wear protective gloves, protective clothing, face protection, eye

P305+P351 - IF IN EYES: Rinse cautiously with water for several minutes.

P337+P313 - If eye irritation persists: Get medical advice/attention.

2.3. Other hazards

No additional information available

SECTION 3: Composition/information on ingredients

3.1. **Substances**

Name : (DL)-Malic acid CAS-No. : 6915-15-7 EC-No. : 230-022-8

Full text of H-statements: see section 16

3.2. **Mixtures**

Not applicable

SECTION 4: First aid measures

4.1. Description of first aid measures

First-aid measures after inhalation : Assure fresh air breathing.

First-aid measures after skin contact : Wash skin with mild soap and water. If skin irritation occurs: Get medical

advice/attention.

: Rinse cautiously with water for several minutes. Remove contact lenses, if First-aid measures after eye contact

present and easy to do. Continue rinsing. Obtain medical attention if pain,

blinking or redness persists.

: Give water to drink. Get medical advice/attention. First-aid measures after ingestion

4.2. Most important symptoms and effects, both acute and delayed

No additional information available

4.3. Indication of any immediate medical attention and special treatment needed

No additional information available

SECTION 5: Firefighting measures

Extinguishing media

Suitable extinguishing media : Dry chemical, CO2, or water spray or regular foam. according to Regulation (EU) 2015/830

5.2. Special hazards arising from the substance or mixture

Hazardous decomposition products in : Under fire conditions, hazardous fumes will be present.

case of fire

5.3. Advice for firefighters

Firefighting instructions : Use water spray or fog for cooling exposed containers. Prevent fire fighting

water from entering the environment.

Protection during firefighting : Do not enter fire area without proper protective equipment, including

respiratory protection.

SECTION 6: Accidental release measures

6.1. Personal precautions, protective equipment and emergency procedures

General measures : Spill should be handled by trained cleaning personnel properly equipped

with respiratory and eye protection.

6.1.1. For non-emergency personnel

Protective equipment : Wear suitable protective clothing.

Measures in case of dust release : Avoid dust formation.

6.1.2. For emergency responders

No additional information available

6.2. Environmental precautions

Prevent product from entering drains. Avoid release to the environment.

6.3. Methods and material for containment and cleaning up

Methods for cleaning up : Sweep up dry powder and dispose properly.

6.4. Reference to other sections

No additional information available

SECTION 7: Handling and storage

7.1. Precautions for safe handling

Precautions for safe handling : Ensure good ventilation of the work station. Handle in accordance with good

industrial hygiene and safety procedures.

7.2. Conditions for safe storage, including any incompatibilities

Storage conditions : Store in a well-ventilated place. Keep container tightly closed. Store at room

temperature.

7.3. Specific end use(s)

For professional use only. Duchefa Biochemie B.V. products are intended only for "in vitro laboratory" research purposes.

SECTION 8: Exposure controls/personal protection

8.1. Control parameters

No additional information available

8.2. Exposure controls

Hand protection:

Safety Data Sheet

according to Regulation (EU) 2015/830

Туре	Material	Permeation	Thickness (mm)	Standard
Gloves	Nitrile rubber (NBR)	6 (> 480 minutes)	0,11	EN 374

Eye protection:

Safety glasses

Skin and body protection:

Wear suitable protective clothing

Respiratory protection:

Filter type P2 (EN 143)

SECTION 9: Physical and chemical properties

9.1. Information on basic physical and chemical properties

Physical state : Solid

Appearance : Powder.

Molecular mass : 134,1 g/mol

Colour : Colourless.

Odour : Weak.

Odour threshold : No data available pH : No data available

pH solution : 2,3 10g/l

Relative evaporation rate : No data available

(butylacetate=1)

Melting point : 128 - 132 °C
Freezing point : No data available
Boiling point : No data available

Flash point : 203 °C

Auto-ignition temperature : No data available
Decomposition temperature : No data available
Flammability (solid, gas) : No data available
Vapour pressure : No data available
Relative vapour density at 20 °C : No data available
Relative density : No data available

Density : 1,6 g/cm³ Solubility : Water: 530 g/l

Log Pow : -1,26

Viscosity, kinematic : No data available
Viscosity, dynamic : No data available
Explosive properties : No data available
Oxidising properties : No data available
Explosive limits : No data available

9.2. Other information

No additional information available

 21/11/2017 Duchefa Biochemie B.V.
 EN (English)
 SDS Ref.: M1315
 4/8

SECTION 10: Stability and reactivity

10.1. Reactivity

No additional information available

10.2. Chemical stability

Stable under normal conditions.

10.3. Possibility of hazardous reactions

No additional information available

10.4. Conditions to avoid

Overheating.

10.5. Incompatible materials

No additional information available

10.6. Hazardous decomposition products

When heated to decomposition, emits dangerous fumes.

SECTION 11: Toxicological information

11.1. Information on toxicological effects

Acute toxicity : Not classified

(DL)-Malic acid (6915-15-7)		
	LD50 oral rat	> 3200 kilogram

Skin corrosion/irritation : Slightly irritating to the skin Serious eye damage/irritation : Causes serious eye irritation.

Respiratory or skin sensitisation : Not classified Germ cell mutagenicity : Not classified Carcinogenicity : Not classified

Reproductive toxicity : Not classified STOT-single exposure : Not classified STOT-repeated exposure : Not classified

Aspiration hazard : Not classified

SECTION 12: Ecological information

12.1. Toxicity

(DL)-Malic acid (6915-15-7)	
LC50 fish 1	295 mg/l Carassius auratus (goldfish)

12.2. Persistence and degradability

(DL)-Malic acid (6915-15-7)	
Biodegradation	> 90 % 3d

12.3. Bioaccumulative potential

(DL)-Malic acid (6915-15-7)	
Log Pow	-1,26

SDS Ref.: M1315

Safety Data Sheet

according to Regulation (EU) 2015/830

12.4. Mobility in soil

No additional information available

12.5. Results of PBT and vPvB assessment

No additional information available

12.6. Other adverse effects

Additional information : Avoid release to the environment. Prevent entry to sewers and public waters

SECTION 13: Disposal considerations

13.1. Waste treatment methods

Waste treatment methods

Dispose in a safe manner in accordance with local/national regulations.
 Dispose of this material and its container at hazardous or special waste collection point.

SECTION 14: Transport information

In accordance with ADR / IATA / IMDG

ADR	IMDG	IATA	
14.1. UN number			
Not applicable	Not applicable	Not applicable	
14.2. UN proper shipping name			
Not applicable	Not applicable	Not applicable	
Not applicable	Not applicable	Not applicable	
14.3. Transport hazard class(es)			
Not applicable	Not applicable	Not applicable	
Not applicable	Not applicable	Not applicable	
14.4. Packing group			
Not applicable	Not applicable	Not applicable	
14.5. Environmental hazards			
Not applicable	Not applicable	Not applicable	
No supplementary information available			

14.6. Special precautions for user

Overland transport

Not applicable

- Transport by sea

Not applicable

- Air transport

Not applicable

14.7. Transport in bulk according to Annex II of Marpol and the IBC Code

Not applicable

SECTION 15: Regulatory information

15.1. Safety, health and environmental regulations/legislation specific for the substance or mixture

15.1.1. EU-Regulations

Safety Data Sheet

according to Regulation (EU) 2015/830

No REACH Annex XVII restrictions

(DL)-Malic acid is not on the REACH Candidate List

(DL)-Malic acid is not on the REACH Annex XIV List

15.1.2. National regulations

Germany

12th Ordinance Implementing the Federal Immission Control Act -

12.BlmSchV

: Is not subject of the 12. BlmSchV (Hazardous Incident Ordinance)

Netherlands

SZW-lijst van kankerverwekkende

stoffen

: The substance is not listed

SZW-lijst van mutagene stoffen NIET-limitatieve lijst van voor de

NIET-limitatieve lijst van voor de voortplanting giftige stoffen –

Borstvoeding

: The substance is not listed: The substance is not listed

NIET-limitatieve lijst van voor de voortplanting giftige stoffen –

Vruchtbaarheid

: The substance is not listed

NIET-limitatieve lijst van voor de voortplanting giftige stoffen –

Ontwikkeling

: The substance is not listed

15.2. Chemical safety assessment

No additional information available

SECTION 16: Other information

Indication of changes:

1.4	Emergency number	Modified	
2	Classification according to Directive 67/548/EEC [DSD] or 1999/45/EC [DPD]	Removed	Obsolete
8.2	Respiratory protection	Modified	Specified P2 protection.
8.2	Hand protection	Modified	Specified material, thickness, et cetera of gloves

Abbreviations and acronyms:

ATE	Acute Toxicity Estimate
ADR	European Agreement concerning the International Carriage of Dangerous Goods by Road
BCF	Bioconcentration factor
CLP	Classification Labelling Packaging Regulation; Regulation (EC) No 1272/2008
DPD	Dangerous Preparations Directive 1999/45/EC
DSD	Dangerous Substances Directive 67/548/EEC
IATA	International Air Transport Association
IMDG	International Maritime Dangerous Goods
LC50	Median lethal concentration
LD50	Median lethal dose
LOAEL	Lowest Observed Adverse Effect Level
NOAEC	No-Observed Adverse Effect Concentration
PBT	Persistent Bioaccumulative Toxic
REACH	Registration, Evaluation, Authorisation and Restriction of Chemicals Regulation (EC) No 1907/2006
SDS	Safety Data Sheet

(DL)-Malic acid Safety Data Sheet

according to Regulation (EU) 2015/830

Data sources

: ECHA (European Chemicals Agency). Supplier. REGULATION (EC) No 1272/2008 OF THE EUROPEAN PARLIAMENT AND OF THE COUNCIL of 16 December 2008 on classification, labelling and packaging of substances and mixtures, amending and repealing Directives 67/548/EEC and 1999/45/EC, and amending Regulation (EC) No 1907/2006.

Full text of H- and EUH-statements:

Eye Irrit. 2	Serious eye damage/eye irritation, Category 2
H319	Causes serious eye irritation.

SDS Biochemicals version 2017

This information is based on our current knowledge and is intended to describe the product for the purposes of health, safety and environmental requirements only. It should not therefore be construed as guaranteeing any specific property of the product

 21/11/2017 Duchefa Biochemie B.V.
 EN (English)
 SDS Ref.: M1315
 8/8