

## FICHE TECHNIQUE

### M1404: MUG TRIHYDRATE

CAS number 6160-80-1

Molecular weight  $C_{16}H_{16}O_9 \cdot 3H_2O = 406.4$

#### Storage condition

- Store dry at 2-8 °C
- Soluble in DMF and DMSO

#### Hazard statements

- **H335** May cause respiratory irritation

#### Precautionary statements

- **P260** Do not breathe dust/fume/gas/mist/vapours/spray
- **P304** IF INHALED:
- **P312** Call a POISON CENTER or doctor/physician if you feel unwell
- **P340** Remove victim to fresh air and keep at rest in a position comfortable for breathing

M1404.0100/ 100 mg

M1404.1000/ 1 g