

Safety Data Sheet

M1415 Version: 2.0

according to Regulation (EU) 2015/830

Date of issue: 05/07/2011 Revision date: 06/11/2017

Supersedes 02/02/2016

SECTION 1: Identification of the substance/mixture and of the company/undertaking

1.1. Product identifier

Product form : Substance

Trade name : MTT (thiazolyl blue)

EC-No. : 206-069-5 CAS-No. : 298-93-1 Product code : M1415

Formula : C18H16N5SBr

Synonyms : Thiazolyl blue tetrazolium bromide / Dimethyl-thiazolyldiphenyl-tetrazolium

bromide / 3-(4,5-Dimethylthiazolyl-2)-2,5-diphenyltetrazolium bromide

Product group : Raw material

1.2. Relevant identified uses of the substance or mixture and uses advised against

1.2.1. Relevant identified uses

Main use category : Professional use

Industrial/Professional use spec : For professional use only. Duchefa Biochemie B.V. products are intended

only for "in vitro laboratory" research purposes.

1.2.2. Uses advised against

No additional information available

1.3. Details of the supplier of the safety data sheet

Supplier

Duchefa Biochemie B.V.
A. Hofmanweg 71
2031 BH Haarlem - The Netherlands
T +31(0)23-5319093 - F +31(0)23-5318027
info@duchefa.nl

1.4. Emergency telephone number

Emergency number : Supplier contact information:

+31(0)23-5319093 (M-F 09:00-17:00) +31(0)6-30109355 (outside office hours)

Organisation/Company	Address	Comment
World Health Organization world	http://apps.who.int/poisoncentres/	Consult website for a local poison centre
directory of poison centres		

SECTION 2: Hazards identification

2.1. Classification of the substance or mixture

Classification according to Regulation (EC) No. 1272/2008 [CLP]

Germ cell mutagenicity, Category 2

H341
Serious eye damage/eye irritation, Category 2

H319
Skin corrosion/irritation, Category 2

H315
Specific target organ toxicity — Single exposure, Category 3,

H335

Respiratory treat irritation

Respiratory tract irritation

Full text of hazard classes and H-statements : see section 16

Safety Data Sheet

according to Regulation (EU) 2015/830

Adverse physicochemical, human health and environmental effects

No additional information available

2.2. Label elements

Labelling according to Regulation (EC) No. 1272/2008 [CLP]

Hazard pictograms (CLP)





GHS07

Signal word (CLP) : Warning

Hazard statements (CLP) : H315 - Causes skin irritation

> H319 - Causes serious eye irritation. H335 - May cause respiratory irritation. H341 - Suspected of causing genetic defects.

: P281 - Use personal protective equipment as required. Precautionary statements (CLP)

P261 - Avoid breathing dust, fume, gas, mist, spray, vapours.

P305+P351 - IF IN EYES: Rinse cautiously with water for several minutes.

2.3. Other hazards

No additional information available

SECTION 3: Composition/information on ingredients

3.1. **Substances**

Name : MTT (thiazolyl blue)

CAS-No. : 298-93-1 EC-No. : 206-069-5

Name	Product identifier	%
MTT (thiazolyl blue)	(CAS-No.) 298-93-1	> 98
	(EC-No.) 206-069-5	

Full text of H-statements: see section 16

3.2. **Mixtures**

Not applicable

SECTION 4: First aid measures

4.1. Description of first aid measures

First-aid measures after inhalation : Move to fresh air. If you feel unwell, seek medical advice.

First-aid measures after skin contact : Seek medical attention if ill effect or irritation develops. Wash off with soap

and plenty of water.

First-aid measures after eye contact : Rinse immediately with plenty of water. Seek medical attention if ill effect or

irritation develops.

First-aid measures after ingestion : Rinse mouth. If you feel unwell, seek medical advice.

4.2. Most important symptoms and effects, both acute and delayed

No additional information available

according to Regulation (EU) 2015/830

4.3. Indication of any immediate medical attention and special treatment needed

No additional information available

SECTION 5: Firefighting measures

5.1. Extinguishing media

Suitable extinguishing media : Alcohol resistant foam. Dry chemical powder. Carbon dioxide (CO2). Water

spray.

Unsuitable extinguishing media : Strong water jet.

5.2. Special hazards arising from the substance or mixture

Hazardous decomposition products in : Under fire conditions, hazardous fumes will be present. COx. NOx. SOx.

case of fire

5.3. Advice for firefighters

Firefighting instructions : Avoid (reject) fire-fighting water to enter environment.

Protection during firefighting : Wear proper protective equipment.

SECTION 6: Accidental release measures

6.1. Personal precautions, protective equipment and emergency procedures

6.1.1. For non-emergency personnel

Protective equipment : Wear suitable protective clothing.

Measures in case of dust release : Use good housekeeping practices to avoid rendering dust airborne.

6.1.2. For emergency responders

No additional information available

6.2. Environmental precautions

Prevent entry to sewers and public waters.

6.3. Methods and material for containment and cleaning up

Methods for cleaning up : Sweep up dry powder and dispose properly.

6.4. Reference to other sections

No additional information available

SECTION 7: Handling and storage

7.1. Precautions for safe handling

Precautions for safe handling : Handle in accordance with good industrial hygiene and safety practice.

7.2. Conditions for safe storage, including any incompatibilities

Storage conditions : Store in dry, cool, well-ventilated area.

Storage temperature : 2 - 8 °C

7.3. Specific end use(s)

For professional use only. Duchefa Biochemie B.V. products are intended only for "in vitro laboratory" research purposes.

06/11/2017 Duchefa Biochemie B.V. EN (English) SDS Ref.: M1415

according to Regulation (EU) 2015/830

SECTION 8: Exposure controls/personal protection

8.1. Control parameters

No additional information available

8.2. Exposure controls

Hand protection:

Туре	Material	Permeation	Thickness (mm)	Standard
Gloves	Nitrile rubber (NBR)	6 (> 480 minutes)	0,11	EN 374

Eye protection:

Safety glasses (to European standard EN 166 or equivalent)

Skin and body protection:

Long sleeved protective clothing. Chemical resistant apron

Respiratory protection:

Wear approved mask. Filter type P3 (EN 143)

SECTION 9: Physical and chemical properties

9.1. Information on basic physical and chemical properties

Physical state : Solid
Appearance : Powder.
Molecular mass : 414,32 g/mol
Colour : Yellow.

Odour : No data available
Odour threshold : No data available
pH : No data available
Relative evaporation rate : No data available

(butylacetate=1)

Melting point : 195 °C

: No data available Freezing point : No data available **Boiling point** Flash point : No data available : No data available Auto-ignition temperature Decomposition temperature : No data available Flammability (solid, gas) : No data available : No data available Vapour pressure Relative vapour density at 20 °C : No data available : No data available Relative density : Water: 10 g/l (20 °C) Solubility

Log Pow : No data available
Viscosity, kinematic : No data available
Viscosity, dynamic : No data available
Explosive properties : No data available
Oxidising properties : No data available
Explosive limits : No data available

06/11/2017 Duchefa Biochemie B.V. EN (English) SDS Ref.: M1415 4/8

Safety Data Sheet

according to Regulation (EU) 2015/830

9.2. Other information

No additional information available

SECTION 10: Stability and reactivity

10.1. Reactivity

Stable under normal conditions of storage, handling and use.

10.2. Chemical stability

Stable under normal conditions.

10.3. Possibility of hazardous reactions

This product is not hazardous.

10.4. Conditions to avoid

Direct sunlight. Heat.

10.5. Incompatible materials

No additional information available

10.6. Hazardous decomposition products

Thermal decomposition generates: COx. NOx. SOx.

SECTION 11: Toxicological information

11.1. Information on toxicological effects

Acute toxicity : Not classified

Skin corrosion/irritation : Causes skin irritation.

Serious eye damage/irritation : Causes serious eye irritation.

Respiratory or skin sensitisation : Not classified

Germ cell mutagenicity : Suspected of causing genetic defects.

Carcinogenicity : Not classified

Reproductive toxicity : Not classified

STOT-single exposure : May cause respiratory irritation.

STOT-repeated exposure : Not classified

Aspiration hazard : Not classified

Other information : See actual entry in RTECS for complete information: XF8060000.

SECTION 12: Ecological information

12.1. Toxicity

No additional information available

12.2. Persistence and degradability

No additional information available

12.3. Bioaccumulative potential

No additional information available

12.4. Mobility in soil

No additional information available

Safety Data Sheet

according to Regulation (EU) 2015/830

12.5. Results of PBT and vPvB assessment

No additional information available

12.6. Other adverse effects

Additional information : Prevent entry to sewers and public waters. Avoid release to the environment

SECTION 13: Disposal considerations

13.1. Waste treatment methods

Waste treatment methods : Dispose in a safe manner in accordance with local/national regulations.

Avoid release to the environment.

SECTION 14: Transport information

In accordance with ADR / IATA / IMDG

ADR	IMDG	IATA		
14.1. UN number				
Not regulated	Not regulated	Not regulated		
14.2. UN proper shipping name	14.2. UN proper shipping name			
Not regulated	Not regulated	Not regulated		
Not regulated	Not regulated	Not regulated		
14.3. Transport hazard class(es)				
Not regulated	Not regulated	Not regulated		
Not regulated	Not regulated	Not regulated		
14.4. Packing group				
Not regulated	Not regulated	Not regulated		
14.5. Environmental hazards				
Not regulated	Not regulated	Not regulated		
No supplementary information available				

14.6. Special precautions for user

. Overland transport

Not regulated

- Transport by sea

Not regulated

- Air transport

Not regulated

14.7. Transport in bulk according to Annex II of Marpol and the IBC Code

Not applicable

SECTION 15: Regulatory information

15.1. Safety, health and environmental regulations/legislation specific for the substance or mixture

15.1.1. EU-Regulations

No REACH Annex XVII restrictions

MTT (thiazolyl blue) is not on the REACH Candidate List

MTT (thiazolyl blue) is not on the REACH Annex XIV List

SDS Ref.: M1415

Safety Data Sheet

according to Regulation (EU) 2015/830

15.1.2. National regulations

Ensure all national/local regulations are observed.

Germany

12th Ordinance Implementing the Federal Immission Control Act -12.BImSchV

: Is not subject of the 12. BlmSchV (Hazardous Incident Ordinance)

Netherlands

SZW-lijst van kankerverwekkende

stoffen

: The substance is not listed

SZW-lijst van mutagene stoffen NIET-limitatieve lijst van voor de voortplanting giftige stoffen -Borstvoeding

: The substance is not listed : The substance is not listed

NIET-limitatieve lijst van voor de voortplanting giftige stoffen -

Vruchtbaarheid

: The substance is not listed

NIET-limitatieve lijst van voor de voortplanting giftige stoffen -

Ontwikkeling

: The substance is not listed

Denmark

Recommendations Danish Regulation : Young people below the age of 18 years are not allowed to use the product Pregnant/breastfeeding women working with the product must not be in direct contact with the product

> The requirements from the Danish Working Environment Authorities regarding work with carcinogens must be followed during use and disposal

Chemical safety assessment 15.2.

No additional information available

SECTION 16: Other information

Indication of changes:

1.4	Emergency number	Modified	
8.2	Hand protection	Modified	Specified material, thickness, et cetera of gloves

Abbreviations and acronyms:

ATE	Acute Toxicity Estimate
ADR	European Agreement concerning the International Carriage of Dangerous Goods by Road
BCF	Bioconcentration factor
CLP	Classification Labelling Packaging Regulation; Regulation (EC) No 1272/2008
DPD	Dangerous Preparations Directive 1999/45/EC
DSD	Dangerous Substances Directive 67/548/EEC
IATA	International Air Transport Association
IMDG	International Maritime Dangerous Goods
LC50	Median lethal concentration
LD50	Median lethal dose
LOAEL	Lowest Observed Adverse Effect Level
NOAEC	No-Observed Adverse Effect Concentration
PBT	Persistent Bioaccumulative Toxic
REACH	Registration, Evaluation, Authorisation and Restriction of Chemicals Regulation (EC) No 1907/2006

Safety Data Sheet

according to Regulation (EU) 2015/830

SDS Safety Data Sheet

Data sources

: ECHA (European Chemicals Agency). Supplier.

Full text of H- and EUH-statements:

Eye Irrit. 2	Serious eye damage/eye irritation, Category 2
Muta. 2	Germ cell mutagenicity, Category 2
Skin Irrit. 2	Skin corrosion/irritation, Category 2
STOT SE 3	Specific target organ toxicity — Single exposure, Category 3, Respiratory tract irritation
H315	Causes skin irritation
H319	Causes serious eye irritation.
H335	May cause respiratory irritation.
H341	Suspected of causing genetic defects.

SDS Biochemicals version 2017

This information is based on our current knowledge and is intended to describe the product for the purposes of health, safety and environmental requirements only. It should not therefore be construed as guaranteeing any specific property of the product

06/11/2017 Duchefa Biochemie B.V. EN (English) SDS Ref.: M1415 8/8