

FICHE TECHNIQUE

N0253: NLN MEDIUM VITAMIN MIXTURE

VITAMINS

	mg/l	µM
D(+)-Biotine	0.05	0.21
Folic Acid	0.50	1.13
L-Glutamine	800.00	5473.83
Gluthatione (reduced)	30.00	97.61
Glycine	2.00	26.64
Myo-Inositol	100.00	554.94
Nicotinic Acid	5.00	40.62
Pyridoxine HCl	0.50	2.43
L-Serine	100.00	951.57
Thiamine HCl	0.50	1.48

Total concentration vitamins: 1038.55 mg/l

N0253.0001/ 1 l
N0253.0010/ 10 l
N0253.0050/ 50 l