

FICHE TECHNIQUE

N0410: NITSCH VITAMIN MIXTURE

Vitamins

mg/l

μM

Biotin

0.05

Folic acid

0.50

Glycine

2.00

myo-Inositol

100.00

Nicotinic acid

5.00

Pyridoxine HCl

0.50

Thiamine HCl

0.50

Nitsch J.P. and Nitsch C., Haploid plants from pollen grains, Science 169, 85 (1969).

Nitsch J.P., Experimental androgenesis in Nicotiana, Phytomorphology 19, 389 (1969).

N0410.0100/ >100 ml

N0410.0250/ >250 ml