

## FICHE TECHNIQUE



### **N0611: NICOTINIC ACID**

CAS number 59-67-6

Molecular weight  $C_6H_5NO_2 = 123.1$

Assay > 99 %

#### **Storage condition**

- Store at room temperature
- Soluble in water (18 g/l / 20 °C)

#### **Hazard statements**

- **H319** Causes serious eye irritation

#### **Precautionary statements**

- **P262** Do not get in eyes
- **P280** Wear protective gloves/protective clothing/eye protection/face protection
- **P305** IF IN EYES:
- **P351** Rinse continuously with water for several minutes

N0611.0100/ 100 g

N0611.0250/ 250 g

N0611.0500/ 500 g

