

Safety Data Sheet

according to Regulation (EU) 2015/830

Date of issue: 15/04/2011 Supersedes 21/11/2017 Revision date 12/01/2018

N0903 Version: 2.1

SECTION 1: Identification of the substance/mixture and of the company/undertaking

1.1. Product identifier

Product form : Substance

Trade name : a-Naphtalene acetic acid

 EC-No.
 : 201-705-8

 CAS-No.
 : 86-87-3

 Product code
 : N0903

 Formula
 : C12H10O2

Synonyms : 2-(1-naphthyl) acetic acid

NAA.

Product group : Raw material

1.2. Relevant identified uses of the substance or mixture and uses advised against

1.2.1. Relevant identified uses

Main use category : Professional use

Industrial/Professional use spec : For professional use only. Duchefa Biochemie B.V. products are intended

only for "in vitro laboratory" research purposes.

1.2.2. Uses advised against

No additional information available

1.3. Details of the supplier of the safety data sheet

Supplier

Duchefa Biochemie B.V.
A. Hofmanweg 71
2031 BH Haarlem - The Netherlands
T +31(0)23-5319093 - F +31(0)23-5318027
info@duchefa.nl

1.4. Emergency telephone number

Emergency number : Supplier contact information:

+31(0)23-5319093 (M-F 09:00-17:00) +31(0)6-30109355 (outside office hours)

Organisation/Company	Address	Comment
World Health Organization world	http://apps.who.int/poisoncentres/	Consult website for a local poison centre
directory of poison centres		

SECTION 2: Hazards identification

2.1. Classification of the substance or mixture

Classification according to Regulation (EC) No. 1272/2008 [CLP]

Acute toxicity (oral), Category 4 H302
Skin corrosion/irritation, Category 2 H315
Serious eye damage/eye irritation, Category 1 H318
Specific target organ toxicity — Single exposure, Category 3, H335

Respiratory tract irritation

Full text of H statements: see section 16

Safety Data Sheet

according to Regulation (EU) 2015/830

Adverse physicochemical, human health and environmental effects

No additional information available

2.2. Label elements

Labelling according to Regulation (EC) No. 1272/2008 [CLP]

Hazard pictograms (CLP) :





GHS07

GHS05

Signal word (CLP) : Danger

Hazard statements (CLP) : H302 - Harmful if swallowed.

H315 - Causes skin irritation

H318 - Causes serious eye damage. H335 - May cause respiratory irritation.

Precautionary statements (CLP) : P261 - Avoid breathing dust, fume, gas, mist, spray, vapours.

P280 - Wear eye protection, face protection, protective clothing, protective

gloves.

P302+P352 - IF ON SKIN: Wash with plenty of soap and water.

P305+P351 - IF IN EYES: Rinse cautiously with water for several minutes. P301+P312+P330 - IF SWALLOWED : Call a POISON CENTER or doctor if

you feel unwell. Rinse mouth.

2.3. Other hazards

No additional information available

SECTION 3: Composition/information on ingredients

3.1. Substances

Name : a-Naphtalene acetic acid

CAS-No. : 86-87-3 EC-No. : 201-705-8

Name	Product identifier	%
a-Naphtalene acetic acid	(CAS-No.) 86-87-3	> 98
	(EC-No.) 201-705-8	

Full text of H-statements: see section 16

3.2. Mixtures

Not applicable

SECTION 4: First aid measures

4.1. Description of first aid measures

First-aid measures after inhalation : Remove person to fresh air and keep comfortable for breathing.

Seek medical attention if ill effect or irritation develops.

First-aid measures after skin contact : Wash skin with plenty of water and soap

Seek medical attention if ill effect or irritation develops.

First-aid measures after eye contact : Rinse immediately with plenty of water

Get immediate medical advice/attention.

12/01/2018 Duchefa Biochemie B.V. EN (English) SDS Ref.: N0903 2/8

Safety Data Sheet

according to Regulation (EU) 2015/830

First-aid measures after ingestion : Rinse mouth thoroughly with water

Obtain emergency medical attention.

4.2. Most important symptoms and effects, both acute and delayed

Symptoms/effects : Diarrhea. Vomiting. Nausea. Irritating.

4.3. Indication of any immediate medical attention and special treatment needed

No additional information available

SECTION 5: Firefighting measures

5.1. Extinguishing media

Suitable extinguishing media : Dry chemical powder.

Carbon dioxide (CO2).

Water spray.

5.2. Special hazards arising from the substance or mixture

Hazardous decomposition products in : Under fire conditions, hazardous fumes will be present

case of fire - COx.

5.3. Advice for firefighters

Firefighting instructions : Use water spray or fog for cooling exposed containers. Avoid (reject) fire-

fighting water to enter environment.

Protection during firefighting : Do not enter fire area without proper protective equipment, including

respiratory protection.

SECTION 6: Accidental release measures

6.1. Personal precautions, protective equipment and emergency procedures

6.1.1. For non-emergency personnel

Protective equipment : Wear suitable protective clothing.

6.1.2. For emergency responders

No additional information available

6.2. Environmental precautions

Prevent entry to sewers and public waters. Avoid release to the environment.

6.3. Methods and material for containment and cleaning up

Methods for cleaning up : Dispose in a safe manner in accordance with local/national regulations.

Sweep up dry powder and dispose properly.

6.4. Reference to other sections

No additional information available

SECTION 7: Handling and storage

7.1. Precautions for safe handling

Precautions for safe handling : Avoid dust formation.

7.2. Conditions for safe storage, including any incompatibilities

Storage conditions : Store in dry, well-ventilated area

Store at room temperature.

 12/01/2018 Duchefa Biochemie B.V.
 EN (English)
 SDS Ref.: N0903
 3/8

Safety Data Sheet

according to Regulation (EU) 2015/830

7.3. Specific end use(s)

For professional use only. Duchefa Biochemie B.V. products are intended only for "in vitro laboratory" research purposes.

SECTION 8: Exposure controls/personal protection

8.1. Control parameters

No additional information available

8.2. Exposure controls

Hand protection:

Туре	Material	Permeation	Thickness (mm)	Standard
Gloves	Nitrile rubber (NBR)	6 (> 480 minutes)	0,11	EN 374

Eye protection:

Safety glasses (to European standard EN 166 or equivalent)

Skin and body protection:

Wear suitable protective clothing

Respiratory protection:

Wear approved mask. Filter type P2 (EN 143)

SECTION 9: Physical and chemical properties

9.1. Information on basic physical and chemical properties

Physical state : Solid

Appearance : Powder.

Molecular mass : 186,21 g/mol

Colour : Off-white
Light yellow.

Odour : Negligible.

Odour threshold : No data available pH : No data available Relative evaporation rate : No data available

(butylacetate=1)

Melting point : 129 - 131.5 °C : No data available Freezing point : No data available **Boiling point** Flash point : No data available Auto-ignition temperature : No data available : No data available Decomposition temperature : No data available Flammability (solid, gas) Vapour pressure : No data available Relative vapour density at 20 °C : No data available : No data available Relative density

Density : 1,1 g/l

Solubility : Soluble in 1N NaOH.

Water: 0,38 g/l

Log Pow : No data available

12/01/2018 Duchefa Biochemie B.V. EN (English) SDS Ref.: N0903 4/8

Safety Data Sheet

according to Regulation (EU) 2015/830

Viscosity, kinematic : No data available
Viscosity, dynamic : No data available
Explosive properties : No data available
Oxidising properties : No data available
Explosive limits : No data available

9.2. Other information

No additional information available

SECTION 10: Stability and reactivity

10.1. Reactivity

Stable under normal conditions.

10.2. Chemical stability

Stable under normal conditions of storage, handling and use.

10.3. Possibility of hazardous reactions

According to process conditions, hazardous decomposition products may be generated.

10.4. Conditions to avoid

- Heat
- Moisture.

10.5. Incompatible materials

- Strong oxidizers

Strong acids.

10.6. Hazardous decomposition products

Thermal decomposition generates:

- COx.

SECTION 11: Toxicological information

11.1. Information on toxicological effects

Acute toxicity : Oral: Harmful if swallowed.

a-Naphtalene acetic acid (86-87-	3)
LD50 oral rat	> 1000 mg/kg
LD50 dermal rabbit	> 5000 mg/kg

Skin corrosion/irritation : Causes skin irritation.

Serious eye damage/irritation : Causes serious eye damage.

Respiratory or skin sensitisation : Not classified Germ cell mutagenicity : Not classified Carcinogenicity : Not classified

Reproductive toxicity : Not classified

STOT-single exposure : May cause respiratory irritation.

STOT-repeated exposure : Not classified

Aspiration hazard : Not classified

Other information : See actual entry in RTECS for complete information: QJ0875000.

12/01/2018 Duchefa Biochemie B.V. EN (English) SDS Ref.: N0903 5/8

Safety Data Sheet

according to Regulation (EU) 2015/830

SECTION 12: Ecological information

12.1. Toxicity

a-Naphtalene acetic acid (86-87-3)	
LC50 fish 1	28 mg/l Oncorhynchus mykiss (Rainbow trout)
EC50 Daphnia 1	180 mg/l EC50 48h - Daphnia magna [mg/l]

12.2. Persistence and degradability

No additional information available

12.3. Bioaccumulative potential

No additional information available

12.4. Mobility in soil

No additional information available

12.5. Results of PBT and vPvB assessment

No additional information available

12.6. Other adverse effects

Additional information : Prevent entry to sewers and public waters. Avoid release to the environment

SECTION 13: Disposal considerations

13.1. Waste treatment methods

Waste treatment methods

: When not empty dispose of this container at hazardous or special waste collection point. Dispose in a safe manner in accordance with local/national regulations. Avoid release to the environment.

SECTION 14: Transport information

In accordance with IMDG / IATA / ADR

ADR	IMDG	IATA	
14.1. UN number			
Not regulated	Not regulated	Not regulated	
14.2. UN proper shipping name			
Not regulated	Not regulated	Not regulated	
Not regulated	Not regulated	Not regulated	
14.3. Transport hazard class(es)			
Not regulated	Not regulated	Not regulated	
Not regulated	Not regulated	Not regulated	
14.4. Packing group			
Not regulated	Not regulated	Not regulated	
14.5. Environmental hazards			
Not regulated	Not regulated	Not regulated	
No supplementary information available			

14.6. Special precautions for user

. Overland transport

Not regulated

Safety Data Sheet

according to Regulation (EU) 2015/830

- Transport by sea

Not regulated

- Air transport

Not regulated

14.7. Transport in bulk according to Annex II of Marpol and the IBC Code

Not applicable

SECTION 15: Regulatory information

15.1. Safety, health and environmental regulations/legislation specific for the substance or mixture 15.1.1. EU-Regulations

No REACH Annex XVII restrictions

a-Naphtalene acetic acid is not on the REACH Candidate List

a-Naphtalene acetic acid is not on the REACH Annex XIV List

15.1.2. National regulations

Ensure all national/local regulations are observed.

Germany

12th Ordinance Implementing the Federal Immission Control Act -

12.BImSchV

: Is not subject of the 12. BlmSchV (Hazardous Incident Ordinance)

Netherlands

SZW-lijst van kankerverwekkende stoffen

SZW-lijst van mutagene stoffen

NIET-limitatieve lijst van voor de voortplanting giftige stoffen – Borstvoeding

NIET-limitatieve lijst van voor de voortplanting giftige stoffen -

Vruchtbaarheid

NIET-limitatieve lijst van voor de voortplanting giftige stoffen – Ontwikkeling

: The substance is not listed

: The substance is not listed: The substance is not listed

: The substance is not listed

: The substance is not listed

Denmark

Recommendations Danish Regulation

: Young people below the age of 18 years are not allowed to use the

product

15.2. Chemical safety assessment

No additional information available

SECTION 16: Other information

Indication of changes:

2.2	Precautionary statements (CLP)	Removed	P304+P340
-----	--------------------------------	---------	-----------

Abbreviations and acronyms:

ATE	Acute Toxicity Estimate
ADR	European Agreement concerning the International Carriage of Dangerous Goods by Road
BCF	Bioconcentration factor
CLP	Classification Labelling Packaging Regulation; Regulation (EC) No 1272/2008
DPD	Dangerous Preparations Directive 1999/45/EC
DSD	Dangerous Substances Directive 67/548/EEC
IATA	International Air Transport Association

Safety Data Sheet

according to Regulation (EU) 2015/830

IMDG	International Maritime Dangerous Goods
LC50	Median lethal concentration
LD50	Median lethal dose
LOAEL	Lowest Observed Adverse Effect Level
NOAEC	No-Observed Adverse Effect Concentration
PBT	Persistent Bioaccumulative Toxic
REACH	Registration, Evaluation, Authorisation and Restriction of Chemicals Regulation (EC) No 1907/2006
SDS	Safety Data Sheet

Data sources

: ECHA (European Chemicals Agency). Supplier.

Full text of H- and EUH-statements:

Acute Tox. 4 (Oral)	Acute toxicity (oral), Category 4
Eye Dam. 1	Serious eye damage/eye irritation, Category 1
Skin Irrit. 2	Skin corrosion/irritation, Category 2
STOT SE 3	Specific target organ toxicity — Single exposure, Category 3, Respiratory tract irritation
H302	Harmful if swallowed.
H315	Causes skin irritation
H318	Causes serious eye damage.
H335	May cause respiratory irritation.

SDS Biochemicals version 2017

This information is based on our current knowledge and is intended to describe the product for the purposes of health, safety and environmental requirements only. It should not therefore be construed as guaranteeing any specific property of the product

12/01/2018 Duchefa Biochemie B.V. EN (English) SDS Ref.: N0903 8/8