

FICHE TECHNIQUE



N0903: A-NAPHTALENE ACETIC ACID

CAS number 86-87-3

Molecular weight $C_{12}H_{10}O_2 = 186.2$

Assay > 98 %

Storage condition

- Store at room temperature
- Soluble in 1 N NaOH

Hazard statements

- **H302** Harmful if swallowed
- **H318** Causes serious eye damage

Precautionary statements

- **P280** Wear protective gloves/protective clothing/eye protection/face protection
- **P301** IF SWALLOWED:
- **P305** IF IN EYES:
- **P312** Call a POISON CENTER or doctor/physician if you feel unwell
- **P330** Rinse mouth
- **P351** Rinse continuously with water for several minutes

N0903.0025/ 25 g

N0903.0050/ 50 g

N0903.0100/ 100 g