

## FICHE TECHNIQUE



### **N0912 : B-NAPHTOXYACETIC ACID**

CAS number 120-23-0

Molecular weight C<sub>12</sub>H<sub>10</sub>O<sub>3</sub> = 202.2

Assay > 97 %

#### **Storage condition**

- Store at room temperature
- Soluble in 1 N NaOH

#### **Hazard statements**

- **H302** Harmful if swallowed
- **H315** Causes skin irritation
- **H319** Causes serious eye irritation
- **H335** May cause respiratory irritation

#### **Precautionary statements**

- **P261** Avoid breathing dust/fume/gas/mist/vapours/spray
- **P280** Wear protective gloves/protective clothing/eye protection/face protection
- **P305** IF IN EYES:
- **P308** IF exposed or concerned:
- **P313** Get medical advice/attention
- **P351** Rinse continuously with water for several minutes

N0912.0025/ 25 g

N0912.0500/ 500 g