

## FICHE TECHNIQUE



### **N1411: NITRO BLUE TETRAZOLIUM (NBT)**

CAS number 298-83-9

Molecular weight  $C_{40}H_{30}Cl_2N_{10}O_6 = 817.6$

Assay > 99 %

Water < 1.5 %

#### **Storage condition**

- Store dry at 2-8 °C
- Soluble in methanol and water
- Protect from light

#### **Hazard statements**

- **H302** Harmful if swallowed
- **H312** Harmful in contact with skin
- **H315** Causes skin irritation
- **H319** Causes serious eye irritation
- **H332** Harmful if inhaled
- **H335** May cause respiratory irritation

#### **Precautionary statements**

- **P261** Avoid breathing dust/fume/gas/mist/vapours/spray
- **P280** Wear protective gloves/protective clothing/eye protection/face protection
- **P301** IF SWALLOWED:
- **P302** IF ON SKIN:
- **P305** IF IN EYES:
- **P312** Call a POISON CENTER or doctor/physician if you feel unwell
- **P351** Rinse continuously with water for several minutes
- **P352** Wash with soap and water

N1411.0100/ 100 mg

N1411.1000/ 1 g